





# QuickSwitch™ Quant Tetramer Kit Starter Guide

For more detailed information, please see the product datasheet.

Product name	Code No.		
	PE-labeled	APC-labeled	BV421-labeled
QuickSwitch™ Quant HLA-A*01:01 Tetramer Kit	TB-7308-K1	TB-7308-K2	TB-7308-K4
QuickSwitch™ Quant HLA-A*02:01 Tetramer Kit	TB-7300-K1	TB-7300-K2	—
QuickSwitch™ Quant HLA-A*03:01 Tetramer Kit	TB-7306-K1	TB-7306-K2	TB-7306-K4
QuickSwitch™ Quant HLA-A*11:01 Tetramer Kit	TB-7304-K1	TB-7304-K2	TB-7304-K4
QuickSwitch™ Quant HLA-A*24:02 Tetramer Kit	TB-7302-K1	—	TB-7302-K4
QuickSwitch™ Quant H-2Kb Tetramer Kit	TB-7400-K1	TB-7400-K2	TB-7400-K4

## Kit components

Component name	Description	Size	Storage
QuickSwitch™ Tetramer	 Exiting Peptide loaded MHC tetramer	500 µL	2-8°C
Peptide Exchange Factor	The proprietary peptide exchange factor	13 µL	-20°C
Magnetic Capture Beads	 Magnetic beads conjugated with a capture antibody specific for tetramers	500 µL	2-8°C
Exiting Peptide Antibody-FITC (25x)	 FITC conjugated antibody reacting against Exiting Peptide	25 µL	2-8°C
Reference Peptide 1 mM	 The peptide for a positive control	13 µL	-20°C
Assay Buffer (10x)	Washing / Dissolving buffer	1.7 mL	2-8°C


## Materials required but not supplied

- Peptides for new specificity tetramers
- Magnetic tray for microplate
- Microtubes
- Micropipettes and disposable tips
- Ultra pure water
- Plate shaker (Labline model 4625 or equivalent)
- Sonicator (Branson Ultrasonic Cleaner Model #B200 or equivalent)
- Flow cytometer
- Vortex
- Round or conical bottom microplates
- Aluminum foil
- DMSO


## WARNINGS AND PRECAUTIONS

1. The Reference Peptide and concentrated Assay Buffer (10x) must be brought to room temperature (20-25°C) before use.
2. When Assay Buffer (10x) is stored at 2-8°C, some reversible precipitation or turbidity may appear. Incubation at 37°C for a few minutes prior to use is recommended to re-solubilize salts.
3. Diluted solutions (antibody and assay buffer) have to be used on the same day as they are prepared. Therefore it is advised to prepare the exact required volumes just before using them.
4. QuickSwitch™ Tetramer and Exiting Peptide Antibody are light sensitive and therefore should be protected from light during storage and during all the steps of the assay.
5. Care should be taken to avoid splashing and well crosscontaminations.
6. All solutions contain sodium azide (≤0.09%) as preservative. If skin or eye contact occurs, wash excessively with water.
7. This current protocol uses a magnet to pellet the beads. It is possible to pellet by centrifugation using a plate holder or by suction using filter plates.

### A Generation of New Specificity Tetramer Using Peptide Exchange



### B Quantification of Peptide Exchange Using Flow Cytometric Sandwich Immunoassay




## A Generation of New Specificity Tetramer Using Peptide Exchange

- 1 Dissolve each lyophilized peptide to be assayed in DMSO to a 10 mM solution.  
(For high affinity peptides, a 1 mM stock solution is a reasonable starting concentration for the assay.)
- 2 Pipet 50  $\mu$ L of QuickSwitch™ Tetramer into a microtube.
- 3 Add 1  $\mu$ L of peptide and mix gently with pipetting.
- 4 Add 1  $\mu$ L of Peptide Exchange Factor and mix gently with pipetting.
- 5 Repeat steps 1 - 4 for each additional peptide.
- 6 Prepare Reference Peptide as follows;
  - (1) Pipet 50  $\mu$ L of QuickSwitch™ Tetramer into a microtube.
  - (2) Add 1  $\mu$ L of Reference Peptide (1 mM) and mix gently with pipetting.
  - (3) Add 1  $\mu$ L of Peptide Exchange Factor and mix gently with pipetting.
- 7 Incubate at least for 4 hours\* at room temperature protected from light.  
(\*Reaction time varies depending on the kit. For details, please refer to the datasheet for each kit.)
- 8 Refrigerate tetramers at 2-8°C protected from light when not used.

## B Quantification of Peptide Exchange Using Flow Cytometric Sandwich Immunoassay

The following assay measures the percentage of original peptide replaced by a competing peptide to help determine whether the resulting tetramer is suitable for antigen-specific CD8<sup>+</sup> T cell staining. All control samples are required for the assay. The QuickSwitch™ Calculator on the website can be downloaded for determining percentages of peptide exchange. [WEB https://www.mblintl.com/quickswitch-peptide-exchange-calculator/](https://www.mblintl.com/quickswitch-peptide-exchange-calculator/)



Control #3 : Beads that have captured QuickSwitch™ Tetramer. FITC conjugated antibody reacts against Exiting Peptide.

Reference Peptide : This peptide typically undergoes about 100% exchange. FITC conjugated antibody does not react against Reference Peptide.

User's Peptide : Peptides selected by users. The figure above shows the flow cytometry data for 50% peptide exchange.

### <Preparation of 1x Assay Buffer>

- 1 Prepare 1x Assay Buffer as follows:  
For 1-5 peptide exchanges, prepare 7.5 mL by mixing 750  $\mu$ L of 10x concentrated Assay Buffer with 6.75 mL of distilled water.  
For 6-8 exchanges, double the volumes.

### <Dispensing Capture Beads>

- Immediately before use, vortex the tetramer capture beads for 30 seconds, followed by a 30-second sonication in a water bath sonicator. If no sonicator is available, vortex an additional 30 seconds.
- Prepare a round or conical-bottom 96 well microtiter plate.
- Pipet 20  $\mu\text{L}$  Magnetic Capture Beads into each of wells #1-#7. Pipet the beads into wells for three essential controls, Reference Peptide, and tetramers generated in Step **A**. The below figure shows the case when three new tetramers are used.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control

### <Tetramer Capture>

- Pipet 5  $\mu\text{L}$  1x Assay Buffer in well #2.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control

- Pipet 5  $\mu\text{L}$  QuickSwitch™ Tetramer in wells #1 and #3.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control

- In wells #4-#7, pipet 5  $\mu\text{L}$  taken from the difference peptide exchange microtubes as described in Step **A**.








Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control

- Shake plate for 45 min. at 550 rpm, protected from light with an opaque cover such as a piece of aluminum foil.



<Rinse>

9 Dispense 150  $\mu$ L of 1x Assay Buffer in wells #1-#7.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							

- 10 Place the plate on a plate magnet for 5 min. to sediment beads.\*
- 11 Flick the plate held tightly to the magnet and blot on a paper towel to minimize cross-contamination of wells.
- 12 After returning the plate upright, vortex for 2 seconds to disperse the beads.








\*It is possible to pellet by centrifugation using a plate holder for 5 min. at 1,000 x g.

<Preparation of 1x Exiting Peptide Antibody>








13 Dilute 25x Exiting Peptide Antibody to 1x as follows: Determine the number (n) of samples to stain with the antibody, including control #2, control #3, Reference Peptides and tetramers generated in Step A. Add one (+1), to account for pipetting errors. In a microtube, pipet (n+1) x 24  $\mu$ L of Assay Buffer and then add (n+1) x 1  $\mu$ L of Exiting Peptide Antibody-FITC. Mix by pipetting.

<Beads incubation with Exiting Peptide Antibody>

14 Pipet 25  $\mu$ L of 1x Exiting Peptide Antibody-FITC prepared in step 13 in wells #2-#7.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							









15 Pipet 25  $\mu$ L of 1x Assay Buffer in well #1.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							









16 Shake plate for 45 min. at 550 rpm, protected from light.

<Rinse>

- 17 Dispense 150  $\mu$ L of 1x Assay Buffer in wells #1-#7.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							









- 18 Place the plate on a plate magnet for 5 min. to sediment beads.\*  
 19 Flick the plate held tightly to the magnet and blot on a paper towel to minimize cross-contamination of wells.  
 20 After returning the plate upright, vortex for 2 seconds to disperse the beads.  
 21 Resuspend beads in 200  $\mu$ L 1x Assay Buffer.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							

\*It is possible to pellet by centrifugation using a plate holder at 1,000 x g for 5 min.

## Preparing Beads Control for flow cytometry assay


- 1 Pipet 200  $\mu$ L 1x Assay Buffer and 5  $\mu$ L Magnetic Capture Beads into well #X.

Well #1	Well #2	Well #3	Well #4	Well #5	Well #6	Well #7	Well #X
Control #1	Control #2	Control #3	Reference Peptide	User's Peptide A	User's Peptide B	User's Peptide C	Beads Control
							

Now they are ready to be analyzed by flow cytometer.


## Flow cytometry analysis

For more detailed information, please see the product datasheet.



## Calculation of results using QuickSwitch™ downloadable calculator

The QuickSwitch™ Calculator can be downloaded for determining percentages of peptide exchange. Please visit the website. [WEB https://www.mblintl.com/quickswitch-peptide-exchange-calculator/](https://www.mblintl.com/quickswitch-peptide-exchange-calculator/)



Percentages of peptide exchange are obtained after MFIs are entered in **the calculation sheet**.



**DOWNLOAD** <https://www.mblintl.com/quickswitch-peptide-exchange-calculator/>

Tetramers bind to T cell receptors (TCR) via three MHC/peptide monomers. Therefore the minimal recommended peptide exchange percentage should be **75%**.

For research use only. Not for use in diagnostic or therapeutic procedures.  
The information is as of June 2022. Please contact us for the latest information. Please read the data sheets carefully before use.  
Copyright © 2022 MEDICAL & BIOLOGICAL LABORATORIES CO., LTD. All Rights Reserved.

2022.06 154043-22061000

**MBL** MEDICAL & BIOLOGICAL  
LABORATORIES CO., LTD.

A JSR Life Sciences Company

SUMITOMO FUDOSAN SHIBADAIMON NICHOME BLDG.  
2-11-8 Shibadaimon, Minato-ku, Tokyo 105-0012 Japan  
TEL: +81-3-6854-3614 E-mail: support@mbi.co.jp  
URL: <https://www.mblbio.com/bio/g/>