

Optimizing MHC Tetramer Staining: Practical Tips - Avoid Common Pitfalls in Tetramer Staining -



Even when following established protocols, a few small adjustments can help improve MHC tetramer staining results.

Key factors such as CD8 antibody selection, the use of appropriate negative controls, and steps to reduce high background signals can all have a significant impact on data quality.

In this issue, we share practical tips based on common questions and troubleshooting cases from researchers working with MHC tetramer reagents. Paying attention to these points can help you obtain clearer and more reliable staining results.

For detailed explanations, recommended conditions, and additional troubleshooting advice, please visit our website.

[Explore MHC Tetramer Optimization and Troubleshooting Tips](#)

Featured Products

Recommended CD8 Clones for MHC Tetramer Staining

For reliable staining with our MHC tetramers, MBL recommends the following anti-CD8 antibody clones.

Certain clones may interfere with tetramer staining, leading to suboptimal results. We strongly recommend clones Hit8a, SFC121Thy2D3 (T8), or MCD8 for human MHC tetramer staining, while for murine experiments CD8 clone selection is even more critical, with clone 53-6.7 to be avoided if possible and clone KT15 strongly recommended.

MBL offers anti-CD8 antibodies validated for compatibility with MHC tetramers:

- [Human CD8: Clone Hit8a](#)
- [Mouse CD8: Clone KT15](#)

Looking for a specific MHC Tetramer?

[→ Search our full range](#)


Need a custom Tetramer?

[→ Contact Us](#)

Or try QuickSwitch™ for rapid in-house preparation

A flexible way to prepare custom tetramers in just a few hours


[→ Learn More](#)



MBL Character, AZ Goes Peptide Mode:Forming a Complex with MHC!

MHC Tetramer

Reagents for direct detection of antigen-specific T cells



Contact us!

Stay tuned for more updates and product information in our upcoming newsletters!
Thank you for choosing MBL for your biological research needs.

For research use only. Not for use in diagnostic or therapeutic procedures. The information is as of June 2026. Please contact us for the latest information. Please read the data sheets carefully before use

Please see the link below for the policy regarding the handling of personal information.
<https://www.mblbio.com/e/privacy/>

MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
SUMITOMO FUDOSAN SHIBADAIMON NICHOME BLDG.
2-11-8 Shibadaimon, Minato-ku, Tokyo 105-0012 Japan
E-mail: inquiry-ruo@mbi.co.jp
<https://www.mblbio.com/e/>

*If you have no recollection of this email, please contact us above.
*Reproducing all or any part of the contents is prohibited.

MEDICAL & BIOLOGICAL LABORATORIES CO., LTD., 2-11-8
Shibadaimon, SUMITOMO FUDOSAN SHIBADAIMON NICHOME
BLDG., Minato-ku, Tokyo, 105-0012, Japan
[Unsubscribe](#) [subscription settings](#)