

# Epitope Tag Antibodies Product Catalog

Wide range of product variation!



---

## Table of Contents

---

Tag Antibodies	P.3
DDDDK-tag	P.4
HA-tag	P.5
His-tag	P.7
Myc-tag	P.9
V5-tag	P.10
mini-AID-tag	P.12
Fluorescent Protein Antibodies	P.13
Other Tag Antibodies	P.18
HRP-Direct	P.22
<i>Smart-IP</i> Tag antibody-conjugated magnetic beads or magnetic agaroses	P.24
Purification Kit & Gel	P.26
FAQs	P.31
Related Antibodies	P.33
Isotype Control Antibodies	
Loading Control Antibodies	
Organelle Marker Antibodies	

### Abbreviation

aff.: affinity purified  
WB: Western blotting, IP: Immunoprecipitation, FCM: Flow cytometry, IC: Immunocytochemistry, IH: Immunohistochemistry,  
IF: Immunofluorescence, RIP: RNP Immunoprecipitation, ChIP: Chromatin Immunoprecipitation,  
Co-IP: Co-Immunoprecipitation, DB: Dot Blotting, NB: Northern Blotting, RNA FISH: RNA Fluorescence *in situ* Hybridization  
\*: Reported in articles, but not tested in MBL. For more details please check the datasheet.

# Japanese High Quality Antibodies

Wide range of product variation,  
from antibodies to purification kits!

## Tag Antibodies

Antigen Target	Antibody	Agarose	Biotin	HRP DirecT Series ▶ P.22 ~	Alexa Fluor® Labeled Antibodies	Magnetic Beads  <i>Smart-IP</i> Series ▶ P.24 ~	Magnetic Agarose  <i>Smart-IP</i> Series ▶ P.24 ~	Purification Kit	
								Agarose ▶ P.26 ~	Magnetic Beads ▶ P.29 ~
<b>DDDDK</b> ▶ P.4	×	×	×	×	×	×	×	×	×
<b>HA</b> ▶ P.5	×	×	×	×	×	×	×	×	×
<b>His</b> ▶ P.7	×	×	×	×	×	×	×	×	
<b>Myc</b> ▶ P.9	×	×	×	×	×	×	×	×	×
<b>V5</b> ▶ P.10	×	×	×	×		×	×	×	×
<b>mini-AID</b> ▶ P.12	×								
<b>GFP</b> ▶ P.13	×	×	×	×	×	×	×		
<b>Renilla GFP</b> ▶ P.14	×								
<b>RFP</b> ▶ P.15	×	×		×		×	×		
<b>S</b> ▶ P.18	×								
<b>E</b> ▶ P.18	×					×	×		
<b>T7</b> ▶ P.18	×	×							
<b>VSV-G</b> ▶ P.18	×	×							
<b>GST</b> ▶ P.18	×			×					
<b>Luciferase</b> ▶ P.19	×								
<b>Renilla Luciferase</b> ▶ P.19	×								

**Go to page 18-21 for more tag antibodies!**

β-galactosidase, MBP, Trx (Thioredoxin), CBD,  
CBP (Calmodulin Binding Protein), Ash-tag, Glu-Glu-tag,  
Strep-tagII, HSV-tag, Digoxigenin (DIG), FITC

HRP-DirecT : HRP is directly labeled to the antibodies for Western blotting.

*Smart-IP* : Antibody-conjugated magnetic beads or magnetic agaroses enabled easy and effective immunoprecipitation of tagged proteins.

Purification kits : Use antibody-conjugated agaroses or magnetic beads for protein purification

Alexa Fluor® is trade mark of Life Technologies. MBL has manufacture and sales the products by receiving the patent license from Life Technologies Corporation in the United States.


## DDDDK-tag

Code No.	Description	Clone	Isotype	Applications	Size
M185-3L	Anti-DDDDK-tag mAb	FLA-1	Mouse IgG2 $\kappa$	WB, IP, FCM, IC	1 mg/1 mL
M185-3LL	Anti-DDDDK-tag mAb	FLA-1	Mouse IgG2 $\kappa$	WB, IP, FCM, IC	1 mg/1 mL x5
M185-3S	Anti-DDDDK-tag mAb	FLA-1	Mouse IgG2 $\kappa$	WB, IP, FCM, IC	50 $\mu$ g/50 $\mu$ L
M185-6	Anti-DDDDK-tag mAb-Biotin	FLA-1	Mouse IgG2 $\kappa$	WB, ELISA	50 $\mu$ L
M185-7	Anti-DDDDK-tag mAb-HRP-Direct <b>HRP-Direct</b>	FLA-1	Mouse IgG2 $\kappa$	WB	200 $\mu$ L
M185-10	Anti-DDDDK-tag mAb-Magnetic Agarose <b>Smart-IP</b>	FLA-1	Mouse IgG2 $\kappa$	IP	100 tests (Gel: 1 mL)
M185-11	Anti-DDDDK-tag mAb-Magnetic Beads <b>Smart-IP</b>	FLA-1	Mouse IgG2 $\kappa$	IP	20 tests (Slurry: 1 mL)
M185-A48	Anti-DDDDK-tag mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	FLA-1	Mouse IgG2 $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L
M185-A59	Anti-DDDDK-tag mAb-Alexa Fluor <sup>®</sup> 594 <b>Alexa Fluor<sup>®</sup></b>	FLA-1	Mouse IgG2 $\kappa$	IC	100 $\mu$ g/100 $\mu$ L
M185-A64	Anti-DDDDK-tag mAb-Alexa Fluor <sup>®</sup> 647 <b>Alexa Fluor<sup>®</sup></b>	FLA-1	Mouse IgG2 $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L


### High sensitivity and high affinity. Can be used in a wide variety of applications.

Applications: WB: 0.1  $\mu$ g/mL (M185-3S, 3L, 3LL)  
 1:2,000 (M185-6)  
 1:2,000-1:5,000 (M185-7)  
 IP: 2  $\mu$ g/sample (M185-3S, 3L, 3LL)  
 10  $\mu$ L/sample (M185-10)  
 50  $\mu$ L/sample (M185-11)  
 FCM: 0.05  $\mu$ g/mL (M185-3S, 3L, 3LL)  
 0.5  $\mu$ g/mL (M185-A48, A64)  
 IC: 0.1  $\mu$ g/mL (M185-3S, 3L, 3LL)  
 0.5-1  $\mu$ g/mL (M185-A48, A59, A64)  
 ELISA: 1:2,000 (M185-6)  
 Specificity: Proteins fused with DDDDK-tag, DYKDDDDK,  
 at N-terminal, Internal, and C-terminal.


#### Western blotting (M185-3L)




#### Western blotting (M185-3L)



#### Immunoprecipitation (M185-3L)



#### Immunocytochemistry (M185-3L)




Code No.	Description	Clone	Isotype	Applications	Size
PM020	Anti-DDDDK-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP, IC	100 $\mu$ L
PM020-7	Anti-DDDDK-tag pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit Ig (aff.)	WB	100 $\mu$ L
PM020-8	Anti-DDDDK-tag pAb-Agarose <b>Agarose</b>	Polyclonal	Rabbit Ig (aff.)	IP	Gel: 200 $\mu$ L


### Reacts with N-terminal, Internal, and C-terminal.

Applications: WB: 1:1,000 (PM020)  
 1:1,000-1:4,000 (PM020-7)  
 IP: 5  $\mu$ L/sample (PM020)  
 20  $\mu$ L/sample (PM020-8)  
 IC: 1:1,000 (PM020)  
 Specificity: Proteins fused with DDDDK-tag, DYKDDDDK,  
 at N-terminal, Internal, and C-terminal.

#### Western blotting (PM020)



#### Immunocytochemistry (PM020)




Code No.	Description	Clone	Isotype	Applications	Size
M180-3	Anti-HA-tag mAb	TANA2	Mouse IgG2b $\kappa$	WB, IP, FCM, IC, RIP*, Co-IP*	200 $\mu$ g/200 $\mu$ L
M180-3S	Anti-HA-tag mAb	TANA2	Mouse IgG2b $\kappa$	WB, IP, FCM, IC	50 $\mu$ g/50 $\mu$ L
M180-6	Anti-HA-tag mAb-Biotin	TANA2	Mouse IgG2b $\kappa$	ELISA	50 $\mu$ L
M180-7	Anti-HA-tag mAb-HRP-Direct <b>HRP-Direct</b>	TANA2	Mouse IgG2b $\kappa$	WB	100 $\mu$ L
M180-10	Anti-HA-tag mAb-Magnetic Agarose <b>Smart-IP</b>	TANA2	Mouse IgG2b $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
M180-11	Anti-HA-tag mAb-Magnetic Beads <b>Smart-IP</b>	TANA2	Mouse IgG2b $\kappa$	IP	20 tests (Slurry: 1 mL)
M180-A48	Anti-HA-tag mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	TANA2	Mouse IgG2b $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L
M180-A59	Anti-HA-tag mAb-Alexa Fluor <sup>®</sup> 594 <b>Alexa Fluor<sup>®</sup></b>	TANA2	Mouse IgG2b $\kappa$	IC	100 $\mu$ g/100 $\mu$ L
M180-A64	Anti-HA-tag mAb-Alexa Fluor <sup>®</sup> 647 <b>Alexa Fluor<sup>®</sup></b>	TANA2	Mouse IgG2b $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L

**High sensitivity. Can be used in a wide variety of applications regardless of the HA-tagged position.**


Applications: WB: 0.1  $\mu$ g/mL (M180-3, 3S)  
 1:5,000-1:10,000 (M180-7)  
 IP: 2  $\mu$ g/sample (M180-3, 3S)  
 10  $\mu$ L/sample (M180-10)  
 50  $\mu$ L/sample (M180-11)  
 FCM: 1  $\mu$ g/mL (M180-3, 3S, A48)  
 2-5  $\mu$ g/mL (M180-A64)  
 IC: 1  $\mu$ g/mL (M180-3, 3S, A48, A59)  
 2-5  $\mu$ g/mL (M180-A64)  
 ELISA: 1:2,000 (M180-6)

Specificity: Proteins fused with HA-tag, YPYDVPDYA,  
 the peptide sequence derived from human influenza hemagglutinin.

**Western blotting (M180-3)**




**Immunoprecipitation (M180-3)**



Lane 1:  
 IP with isotype control (M077-3)  
 Lane 2:  
 IP with Anti-HA-tag mAb (M180-3)

Immunoblotted with 561-7


**Immunocytochemistry (M180-3)**



HA-tagged protein expressed in HeLa

Green: Anti-HA-tag mAb (M180-3)  
 Blue: DAPI


**Flow cytometry (M180-3)**




Open: Anti-HA-tag mAb (M180-3)  
 Closed: Isotype control (M070-3)

**Anti-HA-tag Alexa Fluor<sup>®</sup> 488, 594, 647 (M180-A48, A59, A64)**


Sample: HA-tagged protein expressed in HeLa



Green: Alexa Fluor<sup>®</sup> 488



Red: Alexa Fluor<sup>®</sup> 594



Magenta: Alexa Fluor<sup>®</sup> 647


DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

DDDDK-tag  
HA-tag  
His-tag  
Myc-tag  
V5-tag  
mini-AID-tag  
Fluorescent Protein Antibodies  
Other Tag Antibodies  
HRP-Direct  
Smart-IP  
Purification Kit and Gel  
FAOs  
isotype Control Antibodies  
Loading Control Antibodies  
Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
M132-3	Anti-HA-tag mAb	5D8	Mouse IgG1 $\kappa$	WB, IP	200 $\mu$ g/100 $\mu$ L
M132-10	Anti-HA-tag mAb-Magnetic Agarose <b>Smart-IP</b>	5D8	Mouse IgG1 $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
M132-11	Anti-HA-tag mAb-Magnetic Beads <b>Smart-IP</b>	5D8	Mouse IgG1 $\kappa$	IP	20 tests (Slurry: 1 mL)


Applications: WB: 1  $\mu$ g/mL (M132-3)  
 IP: 1  $\mu$ g/sample (M132-3)  
 10  $\mu$ L/sample (M132-10)  
 50  $\mu$ L/sample (M132-11)

Specificity: Proteins fused with HA-tag, YPYDVPDYA, at N-terminal and C-terminal, including those expressed from pDisplay™ vector.



\*pDisplay™ is a trade mark of Life Technologies Corporation in the United States.

### Western blotting (M132-3)




Code No.	Description	Clone	Isotype	Applications	Size
561	Anti-HA-tag pAb	Polyclonal	Rabbit IgG	WB, IP, IC	100 $\mu$ L
561-5	Anti-HA-tag pAb	Polyclonal	Rabbit IgG	WB, IP, IC	500 $\mu$ L
561-7	Anti-HA-tag pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit IgG	WB	100 $\mu$ L
561-8	Anti-HA-tag pAb-Agarose <b>Agarose</b>	Polyclonal	Rabbit Ig (aff.)	IP	Gel: 200 $\mu$ L
561-A48	Anti-HA-tag pAb-Alexa Fluor® 488 <b>Alexa Fluor®</b>	Polyclonal	Rabbit IgG	FCM, IC	100 $\mu$ L

Applications: WB: 1:1,000 (561, 561-5)  
 1:1,000-1:4,000 (561-7)  
 IP: 1  $\mu$ L/sample (561, 561-5)  
 20  $\mu$ L/sample (561-8)  
 FCM: 1:100 (561-A48)  
 IC: 1:200 (561, 561-5)  
 1:100 (561-A48)


Specificity: Proteins fused with HA-tag, YPYDVPDYA, the peptide sequence derived from human influenza hemagglutinin.

### Western blotting (561)




Sample: HA-tagged protein transfectant

### Immunoprecipitation (561)



Sample: HA-tagged protein transfectant  
 Immunoblotted with 561

### Immunocytochemistry (561)



Anti-HA-tag pAb (561)      Hoechst

HA-Max transfected into BHK cells

Data were provided by Dr. Futoshi Shibasaki, Tokyo Metropolitan Institute of Medical Science.

## His-tag


Code No.	Description	Clone	Isotype	Applications	Size
D291-3	Anti-His-tag mAb	OGHis	Mouse IgG2 $\kappa$	WB, IP, FCM, IC	200 $\mu$ g/200 $\mu$ L
D291-3S	Anti-His-tag mAb	OGHis	Mouse IgG2 $\kappa$	WB, IP, FCM, IC	50 $\mu$ g/50 $\mu$ L
D291-6	Anti-His-tag mAb-Biotin	OGHis	Mouse IgG2 $\kappa$	ELISA	50 $\mu$ g/50 $\mu$ L
D291-7	Anti-His-tag mAb-HRP-Direct <b>HRP-Direct</b>	OGHis	Mouse IgG2 $\kappa$	WB	50 $\mu$ L
D291-8	Anti-His-tag mAb-Agarose <b>Agarose</b>	OGHis	Mouse IgG2 $\kappa$	IP	Gel: 200 $\mu$ L
D291-10	Anti-His-tag mAb-Magnetic Agarose <b>Smart-IP</b>	OGHis	Mouse IgG2 $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
D291-11	Anti-His-tag mAb-Magnetic Beads <b>Smart-IP</b>	OGHis	Mouse IgG2 $\kappa$	IP	20 tests (Slurry: 1 mL)
D291-A48	Anti-His-tag mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	OGHis	Mouse IgG2 $\kappa$	FCM, IC	50 $\mu$ g/50 $\mu$ L
D291-A59	Anti-His-tag mAb-Alexa Fluor <sup>®</sup> 594 <b>Alexa Fluor<sup>®</sup></b>	OGHis	Mouse IgG2 $\kappa$	IC	50 $\mu$ g/50 $\mu$ L
D291-A64	Anti-His-tag mAb-Alexa Fluor <sup>®</sup> 647 <b>Alexa Fluor<sup>®</sup></b>	OGHis	Mouse IgG2 $\kappa$	FCM, IC	50 $\mu$ g/50 $\mu$ L

**High specificity. Can be used for cell staining regardless of the His-tagged position.**

Applications: WB: 0.2  $\mu$ g/mL (D291-3, 3S)  
 1:5,000 (D291-7)  
 IP: 1  $\mu$ g/sample (D291-3, 3S)  
 20  $\mu$ L/sample (D291-8)  
 10  $\mu$ L/sample (D291-10)  
 50  $\mu$ L/sample (D291-11)  
 FCM: 0.2  $\mu$ g/mL (D291-3, 3S)  
 0.5  $\mu$ g/mL (D291-A48, A64)  
 IC: 0.5  $\mu$ g/mL (D291-3, 3S)  
 1  $\mu$ g/mL (D291-A48, A59, A64)  
 ELISA: 0.25-1  $\mu$ g/mL (D291-6)

Specificity: N-terminal, Internal,  
 and C-terminal 6xHis-tagged proteins.


### Western blotting (D291-3)



Lane 1: N-terminal Met-His-tagged protein  
 Lane 2: Internal His-tagged protein  
 Lane 3: Internal His-tagged protein  
 Lane 4: C-terminal His-tagged protein  
 Lane 5: C-terminal His-tagged protein

Met: Methionine


### Immunoprecipitation (D291-3)



Lane 1: IP with Isotype control (M076-3)  
 Lane 2: IP with Anti-His-tag mAb (D291-3)

Immunoblotted with D291-3


### Immunocytochemistry (D291-3)



His-tagged protein expressed in HeLa

Green: Anti-His-tag mAb (D291-3)  
 Blue: DAPI


### Flow cytometry (D291-3)




Open: Isotype control  
 Closed: Anti-His-tag mAb (D291-3)

### Anti-His-tag Alexa Fluor<sup>®</sup> 488, 594, 647 (D291-A48, A59, A64)


Sample: His-tagged protein expressed in HeLa



Green: Alexa Fluor<sup>®</sup> 488  
 Blue: DAPI



Red: Alexa Fluor<sup>®</sup> 594  
 Blue: DAPI



Magenta: Alexa Fluor<sup>®</sup> 647  
 Blue: DAPI

DDDDK-tag  
 HA-tag  
 His-tag  
 Myc-tag  
 V5-tag  
 mini-AID-tag  
 Fluorescent Protein Antibodies  
 Other Tag Antibodies  
 HRP-Direct  
 Smart-IP  
 Purification Kit and Gel  
 FAOs  
 Isotype Control Antibodies  
 Loading Control Antibodies  
 Organelle Marker Antibodies


DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
M136-3	Anti-His-tag mAb	2D8	Mouse IgG2b	IP	100 µg/100 µL

### Great performance in IP.

Applications: IP: 0.5-2.0 µg/sample  
 Specificity: *N*-terminal and *C*-terminal 6xHis-tagged proteins.


#### Immunoprecipitation



Code No.	Description	Clone	Isotype	Applications	Size
M089-3	Anti-His-tag mAb	6C4	Mouse IgG1	WB	100 µg/100 µL

Applications: WB: 1 µg/mL  
 Specificity: *N*-terminal and *C*-terminal 6xHis-tagged proteins.

#### Western blotting




Code No.	Description	Clone	Isotype	Applications	Size
PM032	Anti-His-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 µL
PM032-8	Anti-His-tag pAb-Agarose <span style="background-color: #007bff; color: white; padding: 2px;">Agarose</span>	Polyclonal	Rabbit Ig (aff.)	IP	Gel: 200 µL


### Reacts with *N*-terminal and *C*-terminal.

Applications: WB: 1:1,000 (PM032)  
 IP: 5 µL/sample (PM032)  
 20 µL/sample (PM032-8)  
 Specificity: *N*-terminal and *C*-terminal 6xHis-tagged proteins.

#### Western blotting (PM032)



#### Immunoprecipitation (PM032-8)






Code No.	Description	Clone	Isotype	Applications	Size
M192-3	Anti-Myc-tag mAb	My3	Mouse IgG2b $\kappa$	WB, IP, FCM, IC	200 $\mu$ g/200 $\mu$ L
M192-3S	Anti-Myc-tag mAb	My3	Mouse IgG2b $\kappa$	WB, IP, FCM, IC	50 $\mu$ g/50 $\mu$ L
M192-6	Anti-Myc-tag mAb-Biotin	My3	Mouse IgG2b $\kappa$	WB, ELISA	50 $\mu$ L
M192-7	Anti-Myc-tag mAb-HRP-Direct <b>HRP-Direct</b>	My3	Mouse IgG2b $\kappa$	WB	100 $\mu$ L


**Beats clone 9E10 in sensitivity.**

Applications: WB: 0.1  $\mu$ g/mL (M192-3, 3S)  
 1:10,000 (M192-6, 7)  
 IP: 2  $\mu$ g/sample (M192-3, 3S)  
 FCM: 0.1  $\mu$ g/mL (M192-3, 3S)  
 IC: 0.5  $\mu$ g/mL (M192-3, 3S)  
 ELISA: 1:2,000 (M192-6)  
 Specificity: Proteins fused with Myc-tag, EQKLISEDL

**Western blotting (M192-3)**




Lane 1: Parental cell (293T)  
 Lane 2: N-terminal Myc-tagged protein A/293T  
 Lane 3: C-terminal Myc-tagged protein B/293T  
 Lane 4: C-terminal Myc-tagged protein C/293T  
 Lane 5: C-terminal Myc-tagged protein D/293T  
 Lane 6: Internal Myc-tagged protein E/293T




Sample: C-terminal Myc-tagged protein A/293T

**Immunoprecipitation (M192-3)**




Lane 1: IP with isotype control (M077-3)  
 Lane 2: IP with Anti-Myc-tag mAb (M192-3)  
 Immunoblotted with M047-7


**Flow cytometry (M192-3)**



Open: Anti-Myc-tag mAb (M192-3)  
 Closed: Isotype control (M077-3)



**Immunocytochemistry (M192-3)**




Myc-tagged protein in HeLa  
 Green: Anti-Myc-tag mAb (M192-3)  
 Blue: DAPI

Code No.	Description	Clone	Isotype	Applications	Size
562	Anti-Myc-tag pAb	Polyclonal	Rabbit IgG	WB, IP, IC, IH*, ELISA*, ChIP*	100 $\mu$ L
562-5	Anti-Myc-tag pAb	Polyclonal	Rabbit IgG	WB, IP, IC, IH*, ELISA*, ChIP*	500 $\mu$ L

Applications: WB: 1:1,000  
 IP: 2  $\mu$ L/sample  
 IC: 1:250  
 Specificity: Proteins fused with Myc-tag, EQKLISEDL

**Immunocytochemistry (562)**



Anti-Myc-tag pAb (562)      Hoechst

Myc-tagged Bcl-2 transfected into BHK

Data were provided by Dr. Futoshi Shibasaki, Tokyo Metropolitan Institute of Medical Science.

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAQs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
M047-3	Anti-Myc-tag mAb	PL14	Mouse IgG1 $\kappa$	WB, IP, IC, RIP*	200 $\mu$ g/100 $\mu$ L
M047-6	Anti-Myc-tag mAb-Biotin	PL14	Mouse IgG1 $\kappa$	WB, IC	100 $\mu$ L
M047-7	Anti-Myc-tag mAb-HRP-Direct <b>HRP-Direct</b>	PL14	Mouse IgG1 $\kappa$	WB	100 $\mu$ L
M047-8	Anti-Myc-tag mAb-Agarose <b>Agarose</b>	PL14	Mouse IgG1 $\kappa$	IP	Gel: 200 $\mu$ L
M047-10	Anti-Myc-tag mAb-Magnetic Agarose <b>Smart-IP</b>	PL14	Mouse IgG1 $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
M047-11	Anti-Myc-tag mAb-Magnetic Beads <b>Smart-IP</b>	PL14	Mouse IgG1 $\kappa$	IP	20 tests (Slurry: 1 mL)
M047-A48	Anti-Myc-tag mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	PL14	Mouse IgG1 $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L
M047-A59	Anti-Myc-tag mAb-Alexa Fluor <sup>®</sup> 594 <b>Alexa Fluor<sup>®</sup></b>	PL14	Mouse IgG1 $\kappa$	IC	100 $\mu$ g/100 $\mu$ L
M047-A64	Anti-Myc-tag mAb-Alexa Fluor <sup>®</sup> 647 <b>Alexa Fluor<sup>®</sup></b>	PL14	Mouse IgG1 $\kappa$	FCM, IC	100 $\mu$ g/100 $\mu$ L


### Can be used in wide variety of applications.

Applications: WB: 1  $\mu$ g/mL (M047-3)  
 1:1,000 (M047-6)  
 1:1,000-1:4,000 (M047-7)  
 IP: 5  $\mu$ g/sample (M047-3)  
 20  $\mu$ L/sample (M047-8)  
 10  $\mu$ L/sample (M047-10)  
 50  $\mu$ L/sample (M047-11)

IC: 2  $\mu$ g/mL (M047-3)  
 1:100 (M047-6)  
 5  $\mu$ g/mL (M047-A48, A59)  
 10  $\mu$ g/mL (M047-A64)  
 FCM: 2-5  $\mu$ g/mL (M047-A48)  
 5  $\mu$ g/mL (M047-A64)

Specificity: Proteins fused with Myc-tag, EQKLISEDL.

### Western blotting (M047-3)




### Immunocytochemistry (M047-A48, A59, A64)



Anti-Myc-tag mAb-Alexa Fluor<sup>®</sup>488 (M047-A48)    Anti-Myc-tag mAb-Alexa Fluor<sup>®</sup>594 (M047-A59)    Anti-Myc-tag mAb-Alexa Fluor<sup>®</sup>647 (M047-A64)

Myc-tagged protein expressed in HeLa

### Immunoprecipitation (M047-3)




## V5-tag

Code No.	Description	Clone	Isotype	Applications	Size
M167-3	Anti-V5-tag mAb	1H6	Mouse IgG2 $\alpha$ $\kappa$	WB, IP, IC	100 $\mu$ g/100 $\mu$ L
M167-10	Anti-V5-tag mAb-Magnetic Agarose <b>Smart-IP</b>	1H6	Mouse IgG2 $\alpha$ $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
M167-11	Anti-V5-tag mAb-Magnetic Beads <b>Smart-IP</b>	1H6	Mouse IgG2 $\alpha$ $\kappa$	IP	20 tests (Slurry: 1 mL)

Applications: WB: 1  $\mu$ g/mL (M167-3)  
 IP: 5  $\mu$ g/sample (M167-3)  
 10  $\mu$ L/sample (M167-10)  
 50  $\mu$ L/sample (M167-11)  
 IC: 5  $\mu$ g/mL (M167-3)


Specificity: Proteins fused with V5-tag, GKPIPNPLGLDST, the peptide sequence derived from simian virus 5, at N-terminal and C-terminal.

### Western blotting (M167-3)



Lane1: V5-tagged protein transfectant  
 Lane2: C-terminal V5-tagged protein  
 Lane3: N-terminal V5-tagged protein

### Immunocytochemistry (M167-3)



1: Anti-V5-tag pAb (M167-3)  
 2: GFP own fluorescence  
 3: Merge, 1 and 2  
 4: Transmission light


Code No.	Description	Clone	Isotype	Applications	Size
M215-3	Anti-V5-tag mAb	OZA3	Mouse IgG2bk	WB, IP, FCM, IC	100 µg/100 µL
M215-6	Anti-V5-tag mAb-Biotin	OZA3	Mouse IgG2bk	WB, ELISA	50 µL
M215-7	Anti-V5-tag mAb-HRP-DirectT <b>HRP-DirectT</b>	OZA3	Mouse IgG2bk	WB	100 µL
M215-11	Anti-V5-tag mAb-Magnetic Beads <b>Smart-IP</b>	OZA3	Mouse IgG2bk	IP	20 tests (Slurry: 1 mL)

**Higher performance in WB, IP, and IC as compared to clone 1H6. Can be used in FCM as well!**

Applications: WB: 1 µg/mL (M215-3)  
 1:1,000 (M215-10)  
 IP: 2.5 µg/sample (M215-3)  
 50 µL/sample (M215-11)  
 IC: 1 µg/mL (M215-3)  
 FCM: 0.5 µg/mL (M215-3)


Specificity: Proteins with V5-tag, GKPIPPLLGLDST, the peptide sequence derived from simian virus 5.

### Western blotting (M215-3)




Lane 1: V5-tagged TPO in insect cell culture sup (5 µL/lane)  
 Lane 2: V5-tagged GFP (25 ng/lane)  
 Lane 3: V5-tagged β-galactosidase/HEK293T

### Flow cytometry (M215-3)




Cells: V5-tag GFP/HEK293T  
 Left: Isotype control (M077-3)  
 Right: M215-3

### Immunoprecipitation (M215-3)



Sample: Insect cell culture sup containing V5-tagged TPO  
 Lane 1: Isotype control (M077-3)  
 Lane 2: M215-3  
 Immunoblotted with Anti-V5-tag pAb-HRP-DirectT (PM003-7)

### Immunocytochemistry (M215-3)




Red: M215-3  
 Green: V5-tagged GFP own fluorescence  
 Blue: DAPI  
 Cells: V5-tagged GFP in HeLa transfectant

Code No.	Description	Clone	Isotype	Applications	Size
PM003	Anti-V5-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP, IF*, ChIP*	100 µL
PM003-7	Anti-V5-tag pAb-HRP-DirectT <b>HRP-DirectT</b>	Polyclonal	Rabbit Ig (aff.)	WB	100 µL
PM003-8	Anti-V5-tag pAb-Agarose <b>Agarose</b>	Polyclonal	Rabbit Ig (aff.)	IP	20 tests (Gel: 200 µL)

Applications: WB: 1:2,000 (PM003)  
 1:1,000-1:4,000 (PM003-7)  
 IP: 5 µL/sample (PM003)  
 20 µL/sample (PM003-8)

Specificity: Proteins fused with V5-tag, GKPIPPLLGLDST, the peptide sequence derived from simian virus 5.

### Western blotting / Immunoprecipitation (PM003)



DDDDK-tag  
 HA-tag  
 His-tag  
 Myc-tag  
**V5-tag**  
 mini-AID-tag  
 Fluorescent Protein Antibodies  
 Other Tag Antibodies  
 HRP-DirectT  
 Smart-IP  
 Purification Kit and Gel  
 FAOs  
 Isotype Control Antibodies  
 Loading Control Antibodies  
 Organelle Marker Antibodies

## Auxin-induced rapid protein degradation!


© Both mini-AID-tag and full-length AID-tag can be detected.

Auxin-inducible degron (AID) technology enables to deplete a protein of interest in a half-life less than 30 min by the addition of auxin to culture medium. Auxin includes chemicals such as IAA (indole-3-acetic acid) and NAA (1-naphthalene acetic acid) and includes inhibitors in plants. For the degradation, SCF ubiquitin ligase complex containing the TIR1 protein is activated via the association of auxin with TIR1. Because the core SCF components are conserved in all eukaryotic cells, it is possible to transplant the degradation pathway to yeast and mammalian cells by expressing TIR1. In those cells expressing TIR1, a protein fused with a degron (AID degron) derived from AUX/IAA can be rapidly degraded in the presence of auxin.

mini-AID-tag (8kDa) was developed based on the full-length AID-tag (25kDa). Our original anti-mini-AID-tag, Clone 1E4, can detect mini-AID-tag and full-length AID-tag as well.


\*Expression vector can be obtained by RikenBRC, NBRP-Yeast, or Addgene.

\*The original picture is provided by Dr. Masato Kanemaki from NIG. The picture on the right is edited by MBL.



### ■ Western blotting (M214-3)

Can detect full-length of AID-tag and mini-AID-tag.



### The performance has been proven in articles.

Natsume T *et al.* Rapid Protein Depletion in Human Cells by Auxin-Inducible Degron Tagging with Short Homology Donors. Cell Rep. 15, 210-8 (2016) (PMID:27052166)

#### References

- 1) Nishimura K *et al.* Nature Methods, 6, 917-222(2009) (PMID:19915560)
- 2) Nishimura K, Kanemaki MT. Current Protocols in Cell Biology, 64, 20.9.1-16(2014) (PMID:25181302)

Code No.	Description	Clone	Isotype	Applications	Size
M214-3	Anti-mini-AID-tag mAb	1E4	Mouse IgG2ακ	WB, IP, IC	100 μg/100 μL

## Fluorescent Protein Antibodies

### GFP (Green Fluorescent Protein)


Code No.	Description	Clone	Isotype	Applications	Size
D153-3	Anti-GFP (Green Fluorescent Protein) mAb	RQ2	Rat IgG2aκ	IP, IC	100 µg/100 µL
D153-6	Anti-GFP mAb-Biotin	RQ2	Rat IgG2aκ	ELISA	50 µL
D153-8	Anti-GFP (Green Fluorescent Protein) mAb-Agarose <b>Agarose</b>	RQ2	Rat IgG2aκ	IP, ChIP*, Co-IP*, Purification*	Gel: 200 µL
D153-10	Anti-GFP mAb-Magnetic Agarose <b>Smart-IP</b>	RQ2	Rat IgG2aκ	IP	20 tests (Gel: 200 µL)
D153-11	Anti-GFP (Green Fluorescent Protein) mAb-Magnetic Beads <b>Smart-IP</b>	RQ2	Rat IgG2aκ	IP	20 tests (Slurry: 1 mL)
D153-A48	Anti-GFP (Green Fluorescent Protein) mAb-Alexa Fluor® 488 <b>Alexa Fluor®</b>	RQ2	Rat IgG2aκ	IC	50 µg/50 µL
D153-A59	Anti-GFP (Green Fluorescent Protein) mAb-Alexa Fluor® 594 <b>Alexa Fluor®</b>	RQ2	Rat IgG2aκ	IC	50 µg/50 µL
D153-A64	Anti-GFP (Green Fluorescent Protein) mAb-Alexa Fluor® 647 <b>Alexa Fluor®</b>	RQ2	Rat IgG2aκ	IC	50 µg/50 µL

#### Excellent performance in IP and staining.

Applications: IP: 0.5-2 µg/sample (D153-3)  
 20 µL/sample (D153-8)  
 10 µL/sample (D153-10)  
 50 µL/sample (D153-11)  
 IC: 2 µg/mL (D153-3)  
 2-5 µg/mL (D153-A48, A59, A64)  
 ELISA: 1:5,000/sample (D153-6)


Specificity: Reacts with GFP, EBFP, EGFP, Venus, and Sapphire.

#### Immunoprecipitation (D153-8)




#### Anti-GFP Alexa Fluor® 488, 594, 647 (D153-A48, A59, A64)

Sample: GFP expressed in HeLa



Green: Alexa Fluor® 488,  
 GFP own fluorescence  
 Blue: DAPI



Red: Alexa Fluor® 594  
 Green: DAPI



Magenta: Alexa Fluor® 647  
 Cyan: DAPI


Code No.	Description	Clone	Isotype	Applications	Size
M048-3	Anti-GFP mAb	1E4	Mouse IgG2b	WB, IP, IC, IH	100 µg/100 µL

#### Can be used in a wide variety of applications.

Applications: WB: 1 µg/mL  
 IP: 5 µg/sample  
 IC: 2 µg/mL  
 IH: 10 µg/mL


Specificity: Reacts with GFP, EBFP, SEBFP, ECFP, SECFP, EGFP, SEGFP, cpSEGFP, EYFP, Venus, cpVenus, R-pericam, and Sapphire.

#### Western blotting




Sample: GFP fusion protein transfectant

#### Immunohistochemistry




Anti-GFP mAb (M048-3)



GFP own fluorescence

Paraffin embedded section of GFP transgenic mouse.  
 Magenta: Anti-GFP mAb (M048-3)  
 Green: GFP own fluorescence  
 Cyan and Blue: DAPI

#### Immunohistochemistry



Anti-GFP mAb (M048-3)

Paraffin embedded section of GFP transgenic mouse.

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAQs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
598	Anti-GFP pAb	Polyclonal	Rabbit IgG	WB, IP, IC, IH, ChIP*, Immunoelectron microscopy*	100 µL
598-7	Anti-GFP pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit IgG	WB	100 µL

### The best Anti-GFP.

Applications: WB: 1:1,000-1:5,000 (598)  
1:1,000-1:4,000 (598-7)


IP: 1 µL/sample (598)

IC: 1:500 (598)

IH: 1:1,000-1:2,000 (598)


Specificity: Reacts with GFP, EBFP, SEBFP, ECFP, SECFP, EGFP, SEGFP, cpSEGFP, EYFP, Venus, cpVenus, R-pericam, and Sapphire.

### Immunohistochemistry (598)



GFP Mouse paraffin-embedded section stained with Anti-GFP pAb (598)


### Western blotting (598)



Lane 1: GFP    Lane 4: ECFP\*    Lane 7: SEGFP\*    Lane 10: Venus\*  
Lane 2: EBFP\*    Lane 5: SECFP\*    Lane 8: cpSEGFP\*    Lane 11: cpVenus\*  
Lane 3: SEBFP\*    Lane 6: EGFP\*    Lane 9: EYFP\*    Lane 12: R-pericam\*  
Lane 13: Sapphire\*

\*: Control proteins were provided by RIKEN.

### Immunocytochemistry (598)



Anti-GFP pAb (598)

GFP own fluorescence

GFP expressed in BHK

Data were provided by Dr. Futoshi Shibasaki, Tokyo Metropolitan Institute of Medical Science.

### Reactivity of Anti-GFP and GFP variants

	GFP	EGFP	SEGFP	EBFP	SEBFP	ECFP	SECFP	cpSEGFP	EYFP	Venus	cpVenus	R-pericam	Sapphire
598 (polyclonal)	○	○	○	○	○	○	○	○	○	○	○	○	○
M048-3 (1E4)	○	○	○	○	○	○	○	○	○	○	○	○	○
D153-3 (RQ2)	○	○	N.T.	○	N.T.	○	N.T.	N.T.	N.T.	○	N.T.	N.T.	○


### Renilla GFP

Code No.	Description	Clone	Isotype	Applications	Size
PM073	Anti-Renilla GFP pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP, FCM, IC	100 µL

Applications: WB: 1:1,000  
P: 1 µL/sample  
FCM: 1:2,000  
IC: 1:2,000


Specificity: Reacts with *Renilla* GFP, but not with *Aequorea Victoria* GFP or EGFP.

### Western blotting



Lane 1: *Renilla* GFP/293T  
Lane 2: 293T  
Lane 3: EGFP/293T  
Lane 4: *Aequorea victoria* GFP

### Immunocytochemistry



Magenta: Anti-Renilla GFP pAb (PM073)  
Cyan: DAPI

*Renilla* GFP in HeLa transfectant

## RFP (Red Fluorescent Protein)


Code No.	Description	Clone	Isotype	Applications	Size
M208-3	Anti-RFP mAb Cocktail	1G9, 3G5 (mixed)	Mouse IgG2b <sub>κ</sub> , IgG1 <sub>κ</sub>	WB, IP, FCM, IC	50 μg/50 μL

### © A cocktail of the clone 1G9, good for WB, and clone 3G5, for IC, IP, and FCM.



Applications: WB: 1 μg/mL  
 IP: 5 μg/ sample  
 IC: 1 μg/mL  
 FCM: 0.1-1 μg/mL

Specificity: Reacts with DsRed, mRFP1, mCherry, mOrange, and mPlum.

### ■ Immunocytochemistry (M208-3)



### ■ Western blotting (M208-3)



Code No.	Description	Clone	Isotype	Applications	Size
M204-3	Anti-RFP mAb	1G9	Mouse IgG2b <sub>κ</sub>	WB	100 μg/100 μL
M204-7	Anti-RFP mAb-HRP-Direct	1G9	Mouse IgG2b <sub>κ</sub>	WB	100 μL

### Clone 1G9 is recommended for WB. High sensitivity to detect RFP.


Applications: WB: 1 μg/mL (M204-3)  
 1:5,000 (M204-7)

Specificity: Reacts with DsRed, mRFP1, mCherry, mOrange, and mPlum.

### ■ Western blotting (M204-3)



### ■ Western blotting (M204-3)



### Reactivity of Anti-RFP and RFP variants

	DsRed(tetramer)	mRFP1	mCherry	mOrange	mPlum	mStrawberry	mBanana	mRaspberry
PM005 (polyclonal)	○	○	○	○	○	○	N.T.	N.T.
M155-3 (8D6)	○	○	○	○	○	N.T.	N.T.	N.T.
M165-3 (3G5)	○	○	○	○	○	N.T.	N.T.	N.T.
M204-3 (1G9)	○	○	○	○	○	N.T.	N.T.	N.T.
M208-3 (1G9, 3G5 (mixed))	○	○	○	○	○	N.T.	N.T.	N.T.


DDDDK-tag  
HA-tag  
His-tag  
Myc-tag  
V5-tag  
mini-AID-tag  
Fluorescent Protein Antibodies  
Other Tag Antibodies  
HRP-Direct  
Smart-IP  
Purification Kit and Gel  
FAOs  
Isotype Control Antibodies  
Loading Control Antibodies  
Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
M165-3	Anti-RFP mAb	3G5	Mouse IgG1 $\kappa$	IP, FCM, IC	100 $\mu$ g/100 $\mu$ L
M165-8	Anti-RFP mAb-Agarose <b>Agarose</b>	3G5	Mouse IgG1 $\kappa$	IP	Gel: 200 $\mu$ L
M165-11	Anti-RFP mAb-Magnetic Beads <b>Smart-IP</b>	3G5	Mouse IgG1 $\kappa$	IP	20 tests (Slurry: 1 mL)

### Excellent performance in IP, even the single RFP variant.

Applications: IP: 5  $\mu$ g/sample (M165-3)  
 20  $\mu$ L/sample (M165-8)  
 10  $\mu$ L/sample (M165-10)  
 50  $\mu$ L/sample (M165-11)  
 IC: 1  $\mu$ g/mL (M165-3)  
 FCM: 0.1-1  $\mu$ g/mL (M165-3)  
 Specificity: Reacts with DsRed, mRFP1, mCherry, mOrange, and mPlum.


#### Immunoprecipitation (M165-8)




Code No.	Description	Clone	Isotype	Applications	Size
M155-3	Anti-RFP mAb	8D6	Mouse IgG1 $\kappa$	WB, IC	100 $\mu$ g/100 $\mu$ L

Applications: WB: 1  $\mu$ g/mL  
 IC: 10  $\mu$ g/mL  
 Specificity: Reacts with DsRed, mRFP1, mCherry, mOrange, and mPlum.

#### Western blotting




#### Immunocytochemistry



Code No.	Description	Clone	Isotype	Applications	Size
PM005	Anti-RFP pAb	Polyclonal	Rabbit IgG	WB, IC, IH*	100 $\mu$ L
PM005-7	Anti-RFP pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit IgG	WB	100 $\mu$ L

Applications: WB: 1:1,000 (PM005)  
 1:1,000-1:4,000 (PM005-7)  
 IC: 1:500 (PM005)  
 Specificity: Reacts with DsRed, mRFP, mCherry, mOrange, mPlum, mStrawberry, and tdTomato.

#### Western blotting (PM005)






## Amalgam Fluorescent Protein Antibodies

### © Original Fluorescent Proteins

Code No.	Description	Clone	Isotype	Applications	Size	Cross-reactivity by WB
M102-3M	Anti-monomeric Azami-Green1 mAb	2F11	Mouse IgG1 $\kappa$	WB	100 $\mu$ g/100 $\mu$ L	mAG1
PM052M	Anti-monomeric Azami-Green1 pAb	Polyclonal	Rabbit Ig(aff.)	WB, IP, IC, IH	100 $\mu$ L	mAG1
PM011M	Anti-Azami-Green pAb	Polyclonal	Rabbit IgG	WB	100 $\mu$ L	AG, mAG1
M118-3M	Anti-Dronpa-Green mAb	2F6	Mouse IgG2b	IP	100 $\mu$ g/100 $\mu$ L	
M106-3M	Anti-Kaede mAb	2F4	Mouse IgG1 $\kappa$	IP	100 $\mu$ g/100 $\mu$ L	
M125-3M	Anti-Kaede mAb	3B1	Mouse IgG1	WB	100 $\mu$ g/100 $\mu$ L	Kaede
PM012M	Anti-Kaede pAb	Polyclonal	Rabbit Ig(aff.)	WB	100 $\mu$ L	Kaede
M126-3M	Anti-monomeric Keima-Red mAb	2F7	Mouse IgG2a	WB	100 $\mu$ g/100 $\mu$ L	mKeima-Red
M182-3M	Anti-Keima-Red mAb	1C3	Mouse IgG1 $\kappa$	WB	100 $\mu$ g/100 $\mu$ L	mKeima-Red, dKeima-Red, dKeima570
M128-3M	Anti-Kikume Green-Red mAb	5B3	Mouse IgG2b	WB	100 $\mu$ g/100 $\mu$ L	KikGR, mKikGR
M104-3M	Anti-monomeric Kusabira-Orange1 mAb	1H7	Mouse IgG1 $\kappa$	WB	100 $\mu$ g/100 $\mu$ L	mKO1, mKO2, mKG, mKG-O, mKOkappa
M168-3M	Anti-monomeric Kusabira-Orange2 mAb	3B3	Mouse IgG1 $\kappa$	WB, IP, IC, IH	100 $\mu$ g/100 $\mu$ L	mKO2, mKG, mKG-O, mKOkappa
PM051M	Anti-monomeric Kusabira-Orange2 pAb	Polyclonal	Rabbit Ig(aff.)	WB, IP, IC, IH	100 $\mu$ L	KO1, mKO1, mKO2, mKG, mKG-O, mKOkappa
M116-3M	Anti-Midoriishi-Cyan mAb	2C1	Mouse IgG2b	IP	100 $\mu$ g/100 $\mu$ L	
M130-3M	Anti-Midoriishi-Cyan mAb	5B7	Mouse IgG1	WB	100 $\mu$ g/100 $\mu$ L	MiCy, mMicy
M148-3M	Anti-monomeric Kusabira-Green <i>N</i> -terminal Fragment mAb	1E6	Mouse IgG2b	WB	100 $\mu$ g/100 $\mu$ L	mKO1, mKO2, mKG
M149-3M	Anti-monomeric Kusabira-Green <i>C</i> -terminal Fragment mAb	21B10	Mouse IgG2a	WB	100 $\mu$ g/100 $\mu$ L	mKO2, mKG


#### ■ Immunohistochemistry (PM052M)



Immunohistochemical detection of mAG1 on frozen section of B6. Cg-Tg (Fucci) 504Bsi mouse embryonic brain (E13) with PM052M (1) and Fucci-S/G<sub>2</sub>/M Green own fluorescence (2).

Cg-Tg (Fucci) 596Bsi mouse was provided by RIKEN.

#### ■ Immunocytochemistry (PM051M)



1: Anti-mKO2  
2: Fucci-G1 Orange  
3: Transmission light

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAQs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies


## S-tag

Code No.	Description	Clone	Isotype	Applications	Size
PM021	Anti-S-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 $\mu$ L

Applications: WB: 1:1,000 (PM021)  
IP: 5  $\mu$ L/sample (PM021)  
20  $\mu$ L/sample (PM021-8)

Specificity: Proteins fused with S-tag, KETAAAKFERQHIDS, the peptide sequence derived from Pancreatic RNase A.

### Western blotting (PM021) Immunoprecipitation (PM021)




## E-tag

Code No.	Description	Clone	Isotype	Applications	Size
PM070	Anti-E-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 $\mu$ L

Applications: WB: 1:1,000  
IP: 2  $\mu$ L/sample

Specificity: Proteins fused with E-tag, GAPVYPDPLEPR.

### Western blotting Immunoprecipitation




Code No.	Description	Clone	Isotype	Applications	Size
M198-9	Anti-E-tag mAb-Magnetic Beads <b>Smart-IP</b>	21D11	Mouse IgG2a $\kappa$	IP	20 tests (Slurry: 1 mL)
M198-10	Anti-E-tag mAb-Magnetic Agarose <b>Smart-IP</b>	21D11	Mouse IgG2a $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)

Applications: IP: 50  $\mu$ L/sample (M198-9)  
10  $\mu$ L/sample (M198-10)

Specificity: Proteins fused with E-tag, GAPVYPDPLEPR.

### Immunoprecipitation (M198-10)




## T7-tag

Code No.	Description	Clone	Isotype	Applications	Size
PM022	Anti-T7-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP, ChIP*	100 $\mu$ L
PM022-8	Anti-T7-tag pAb-Agarose <b>Agarose</b>	Polyclonal	Rabbit Ig (aff.)	IP, Purification, RNA pull-down*	Gel: 200 $\mu$ L

Applications: WB: 1:1,000 (PM022)  
IP: 5  $\mu$ L/sample (PM022)  
20  $\mu$ L/sample (PM022-8)  
Specificity: Proteins fused with T7-tag, MASMTGGQQMG, the peptide sequence derived from a capsid protein of T7 bacteriophage.

### Western blotting (PM022) Immunoprecipitation (PM022)




## VSV-G-tag

Code No.	Description	Clone	Isotype	Applications	Size
563	Anti-VSV-G-tag pAb	Polyclonal	Rabbit IgG	WB, IC, IH*	100 $\mu$ L
563-8	Anti-VSV-G-tag pAb-Agarose <b>Agarose</b>	Polyclonal	Rabbit Ig (aff.)	IP	Gel: 200 $\mu$ L

Applications: WB: 1:100 (563)  
IP: 20  $\mu$ L/sample (563-8)  
IC: 1:100 (563)  
Specificity: Reacts with VSV-G-tag, YTDIEMNRLGK.

### Western blotting (563)




## GST (Glutathione-S-transferase)

Code No.	Description	Clone	Isotype	Applications	Size
PM013	Anti-GST-tag pAb	Polyclonal	Rabbit IgG	WB, IP	100 $\mu$ L
PM013-7	Anti-GST-tag pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit Ig (aff.)	WB	50 $\mu$ L

Applications: WB: 1:1,000 (PM013)  
1:5,000 (PM013-7)  
IP: 5  $\mu$ L/sample (PM013)  
Specificity: Proteins fused with GST-tag.

### Western blotting (PM013)




Code No.	Description	Clone	Isotype	Applications	Size
M209-3	Anti-GST-tag mAb	GT5	Mouse IgG1 $\kappa$	WB, IP	100 $\mu$ g/ 100 $\mu$ L

### Clone GT5 is recommended for Western blotting.

Applications: WB: 1  $\mu$ g/mL  
IP: 2.5  $\mu$ g/sample  
Specificity: Proteins fused with GST-tag.


### Western blotting (M209-3)



Code No.	Description	Clone	Isotype	Applications	Size
M071-3	Anti-GST-tag mAb	3B2	Mouse IgG2bx	WB	100 µg/ 100 µL

Application: WB: 1 µg/mL  
Specificity: Proteins fused with GST-tag.

#### Western blotting




## Luciferase

Code No.	Description	Clone	Isotype	Applications	Size
PM016	Anti-Luciferase pAb	Polyclonal	Rabbit IgG	WB, IP, IC, IH	100 µL


Applications: WB: 1:1,000  
IP: 5 µL/sample  
IC: 1:200  
IH: 1:100

Specificity: Reacts with firefly luciferase but not with *Renilla* luciferase.

#### Western blotting



#### Immunoprecipitation




Code No.	Description	Clone	Isotype	Applications	Size
M095-3	Anti-Luciferase mAb	2D4	Mouse IgG1	WB, IP	100 µg/ 100 µL


Applications: WB: 1 µg/mL  
IP: 5 µg/sample

Specificity: Reacts with firefly luciferase but not with *Renilla* luciferase.

#### Western blotting



#### Immunoprecipitation




## Renilla Luciferase

Code No.	Description	Clone	Isotype	Applications	Size
PM047	Anti-Renilla Luciferase pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP, IC	100 µL

Applications: WB: 1:1,000  
IP: 2 µL/sample  
IC: 1:100-1:200  
IC (Paraffin embedded sections): 1:1,000  
(Heat treatment is necessary.)

Specificity: Reacts with *Renilla* luciferase but not with firefly luciferase.

#### Western blotting




## β-galactosidase

Code No.	Description	Clone	Isotype	Applications	Size
PM049	Anti-β-galactosidase pAb	Polyclonal	Rabbit IgG	WB, IP, IC, IH	100 µL

Applications: WB: 1:1,000  
IP: 1 µL/sample  
IC: 1:100  
IH: 1:200

Specificity: β-Galactosidase fusion proteins.

#### Immunohistochemistry (frozen section)



Frozen section of Lewis lung carcinoma xenograft



This data was kindly provided by Dr. Takashi Minami. (Laboratory for Systems Biology and Medicine at RCAST, University of Tokyo)

Code No.	Description	Clone	Isotype	Applications	Size
M203-3	Anti-β-galactosidase mAb	6F4	Rat IgG2ak	WB, IP, IC, IH	100 µg/ 100 µL

Applications: WB: 0.1 µg/mL  
IP: 2 µg/sample  
IC: 2 µg/mL  
IH: 5 µg/mL  
(frozen section)

Specificity: β-Galactosidase fusion proteins.

#### Immunohistochemistry (frozen section)



Tissue: Mouse skeletal muscle (β-galactosidase knock-in)

Green: Anti-β-galactosidase mAb (M203-3)  
Red: Anti-mEndoglin pAb (R&D Systems)  
Blue: DAPI


Data were kindly provided by Ms. Chihito Makihara and Dr. Takashi Minami. (Division of Vascular Biology, RCAST, The University of Tokyo)

Code No.	Description	Clone	Isotype	Applications	Size
M094-3	Anti-β-galactosidase mAb	5A3	Mouse IgG1	WB, IP, FCM*, IC, IH	100 µg/ 100 µL

Applications: WB: 1 µg/mL  
IP: 1 µg/sample  
IC: 5 µg/mL  
IH: 10 µg/mL

Specificity: β-Galactosidase fusion proteins.

#### Immunocytochemistry




Immunoperoxidase staining of β-galactosidase expressed in 293T

## MBP (Maltose Binding Protein)


Code No.	Description	Clone	Isotype	Applications	Size
M091-3	Anti-MBP (Maltose Binding Protein) mAb	1G12	Mouse IgG3	WB, IP	100 µg/ 100 µL

Applications: WB: 1 µg/mL  
IP: 2 µg/sample  
Specificity: Proteins fused with MBP.

#### Western blotting



#### Immunoprecipitation



DDDDK-tag  
HA-tag  
His-tag  
Myc-tag  
V5-tag  
mini-AID-tag  
Fluorescent Protein Antibodies  
Other Tag Antibodies  
HRP-Direct  
Smart-IP  
Purification Kit and Gel  
FAOs  
Isotype Control Antibodies  
Loading Control Antibodies  
Organelle Marker Antibodies


DDDDK-tag  
HA-tag  
His-tag  
Myc-tag  
V5-tag  
mini-AID-tag  
Fluorescent Protein Antibodies  
Other Tag Antibodies  
HRP-Direct  
Smart-IP  
Purification Kit and Gel  
FAOs  
Isotype Control Antibodies  
Loading Control Antibodies  
Organelle Marker Antibodies

## Trx (Thioredoxin)

Code No.	Description	Clone	Isotype	Applications	Size
M013-3	Anti-Thioredoxin (Trx-tag) mAb	2C9	Mouse IgG1κ	WB	100 µg/ 100 µL

Application: WB: 1 µg/mL  
Specificity: Trx and Trx-tag fusion proteins.

### Western blotting




Sample: Thioredoxin fusion protein expressed cell lysate

## CBD (Chitin Binding Domain)

Code No.	Description	Clone	Isotype	Applications	Size
PM015	Anti-CBD (Chitin Binding Domain) pAb	Polyclonal	Rabbit IgG	WB	100 µL

Application: WB: 1:1,000  
Specificity: Proteins fused with CBD-tag.  
TTNPGVSAWQVNTAYTAGQLVIYNGKTYK.

### Western blotting




## CBP (Calmodulin Binding Protein)


Code No.	Description	Clone	Isotype	Applications	Size
PM071	Anti-Calmodulin Binding Protein-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 µL

Application: WB: 1:1,000  
IP: 2-5 µL/sample  
Specificity: Proteins fused with CBD-tag.  
TTNPGVSAWQVNTAYTAGQLVIYNGKTYK.

### Western blotting



### Immunoprecipitation



Immunoblotted with PM071


## Ash-tag

Code No.	Description	Clone	Isotype	Applications	Size
M223-3	Anti-Ash-tag mAb	FLP1C15-2	Mouse IgG1	WB	100 µg/ 100 µL

Application:

WB: 5 µg/mL  
\*Reacts with the code AM8011M and code AM8012M but not with code SI-8010, SI-8011, SI-8020, and SI-8021.

### Western blotting



## Glu-Glu-tag


Code No.	Description	Clone	Isotype	Applications	Size
BRP0001	Anti-Glu-Glu-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 µL

Applications:

WB: 1:1,000  
IP: 5 µL/sample

Specificity: Proteins fused with Glu-Glu-tag, EEEEEYMPME.

### Immunoprecipitation



## Strep-tagII


Code No.	Description	Clone	Isotype	Applications	Size
M211-3	Anti-Strep-tag II mAb	4F1	Mouse IgG2aκ	WB, IP, IC	100 µg/ 100 µL

Applications:

WB: 1 µg/mL  
IP: 2 µg/sample  
IC: 1 µg/mL

Specificity: Proteins fused with Strep-tag II, WSHPPQFEK.

### Immunocytochemistry



## HSV-tag


Code No.	Description	Clone	Isotype	Applications	Size
BRP0002	Anti-HSV-tag pAb	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 µL

Applications:

WB: 1:1,000  
IP: 5 µL/sample

Specificity: Proteins fused with HSV-tag, QPELAPDPED.

### Immunoprecipitation



## Digoxigenin (DIG) , FITC


### Antibodies for labeled nucleic acid probe detection.

Digoxigenin (DIG) and fluorescein like FITC are haptens that can be easily conjugated or integrated into biomolecules such as oligonucleotides and proteins.

These haptens are detected by specific antibodies in wide variety of applications. MBL's monoclonal antibodies against DIG and FITC are proprietary our original clones with high affinity and specificity. They can be used for non-radioactive immunoassays, Northern/Southern blot analysis, and *in situ* hybridization to detect DIG or FITC integrated DNA/RNA probes.


### Anti-Digoxigenin (DIG) mAb (M227-3)

#### Northern blotting



Sample: Total RNA extracted from HEK293T cells  
 Probe: DIG-labeled MALAT-1 ncRNA (RefSeq ID: NR\_002819.3, region 6641-7113)  
 Antibody: M227-3, 1 µg/mL  
 Reagent: DIG Wash and Block Buffer Set (Sigma-Aldrich; code no. 11585762001)  
 Lane 1: 500 ng of total RNA  
 Lane 2: 100 ng of total RNA  
 Lane 3: 20 ng of total RNA


#### RNA Fluorescence *in situ* Hybridization



Cells: HeLa, Probe: MALAT-1 ncRNA (RefSeq ID: NR\_002819.3, region 6641-7113)  
 Antibody: M227-3, 1 µg/mL


### Anti-FITC mAb (M228-3)

#### Dot blotting



Sample: FITC-labeled RNA synthesized by *in vitro* transcription from *lacZ*-encoding cDNA (RefSeq ID: NC\_007779.1, region 363130-364149)  
 <Immunoblot>  
 Lane 1: M228-3  
 Lane 2: Isotype control (M075-3)

#### RNA Fluorescence *in situ* Hybridization



Cells: HeLa,  
 Probe: FITC-labeled MALAT-1 ncRNA (RefSeq ID: NR\_002819.3, region 6641-7113)  
 Antibody: M228-3, 0.5 µg/mL  
 Red: MALAT-1, Blue: DAPI

Code No.	Description	Clone	Isotype	Applications	Size
M227-3	Anti-Digoxigenin (DIG) mAb	8-10	Mouse IgG1κ	WB, IP, ELISA, DB, NB, RNA FISH	100 µg/100 µL
M228-3	Anti-FITC mAb	47-11	Mouse IgG1κ	WB, IP, ELISA, DB, NB, RNA FISH	100 µg/100 µL

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

## HRP-Direct




- Save time and get clear data with MBL's HRP-Direct series.
- No more secondary antibodies for Western blotting!

"HRP-Direct series" is a series of epitope tag antibodies labeled directly with HRP.


MBL has succeeded to develop the HRP labeled antibodies with higher sensitivity and lower background.

**Quick** 30 min for Western blotting.

**Conventional method** 2.5 hrs



**High sensitivity** Can get high sensitivity and save antibody!




**Clear** No cross-reactivity with antibodies used for immunoprecipitation eliminating their heavy and light chain bands in the blots.

Recommended especially when the band size of the target protein overlaps with either heavy chain or light chain.


### Immunoprecipitation (IP)



### Western blotting




■ Anti-DDDDK-tag pAb-HRP-Direct (PM020-7)




IP Antibodies:  
Lane 1: Normal Rabbit IgG-Agarose  
Lane 2: Anti-DDDDK-tag-Agarose

■ Anti-Myc-tag mAb-HRP-Direct (M047-7)




IP Antibodies:  
Lane 1: Isotype control-Agarose  
Lane 2: Anti-Myc-tag-Agarose

■ Anti-HA-tag pAb-HRP-Direct (561-7)




IP Antibodies:  
Lane 1: Normal Rabbit IgG-Agarose  
Lane 2: Anti-HA-tag-Agarose

■ Anti-V5-tag pAb-HRP-Direct (PM003-7)




IP Antibodies:  
Lane 1: Normal Rabbit IgG-Agarose  
Lane 2: Anti-V5-tag-Agarose

■ Anti-His-tag mAb-HRP-Direct (D291-7)



IP Antibody: Anti-His-tag-Agarose  
Lane 1: Conventional method  
Lane 2: Anti-His-tag mAb-HRP-Direct

■ Anti-GFP pAb-HRP-Direct (598-7)  
Anti-RFP pAb-HRP-Direct (PM005-7)



Sample:  
Lane 1, 3: mRFP in 293T  
Lane 2, 4: GFP in 293T  
WB:  
Lane 1, 2: Anti-RFP pAb-HRP-Direct  
Lane 3, 4: Anti-GFP pAb-HRP-Direct

HRP-Direct series

Code No.	Description	Clone	Isotype	Applications	Size
PM020-7	Anti-DDDDK-tag pAb-HRP-Direct	Polyclonal	Rabbit Ig (aff.)	WB	100 µL
M185-7	Anti-DDDDK-tag mAb-HRP-Direct	FLA-1	Mouse IgG2a $\kappa$	WB	200 µL
561-7	Anti-HA-tag pAb-HRP-Direct	Polyclonal	Rabbit IgG	WB	100 µL
M180-7	Anti-HA-tag mAb-HRP-Direct	TANA2	Mouse IgG2b $\kappa$	WB	100 µL
D291-7	Anti-His-tag mAb-HRP-Direct	OGHis	Mouse IgG2a $\kappa$	WB	50 µL
M192-7	Anti-Myc-tag mAb-HRP-Direct	My3	Mouse IgG2b $\kappa$	WB	100 µL
M047-7	Anti-Myc-tag mAb-HRP-Direct	PL14	Mouse IgG1 $\kappa$	WB	100 µL
PM003-7	Anti-V5-tag pAb-HRP-Direct	Polyclonal	Rabbit Ig (aff.)	WB	100 µL
M215-7	Anti-V5-tag mAb-HRP-Direct	OZA3	Mouse IgG2b $\kappa$	WB	100 µL
598-7	Anti-GFP pAb-HRP-Direct	Polyclonal	Rabbit IgG	WB	100 µL
M204-7	Anti-RFP mAb-HRP-Direct	1G9	Mouse IgG2b $\kappa$	WB	100 µL
PM005-7	Anti-RFP pAb-HRP-Direct	Polyclonal	Rabbit IgG	WB	100 µL
PM013-7	Anti-GST-tag pAb-HRP-Direct	Polyclonal	Rabbit Ig (aff.)	WB	50 µL
8480	Loading Control Set ( $\beta$ -Actin, $\alpha$ -Tubulin, GAPDH) HRP-Direct	—	—	—	3 vials 50 µL each

## Tag antibody-conjugated magnetic beads or magnetic agarose

© **Smart-IP** is a series of tag-antibodies conjugated with magnetic beads or magnetic agarose.

**Quick**

Because of using magnetic separation, no need of centrifugal separation.

**High yield**

Magnetic beads are tightly collected on the side of tubes so that washing does not affect the yield quantum.




\*Sold separately Code No. 3190

Type of beads	Applications	Size of beads	Schematic diagram of beads	Product images	Product features	Recommended to who...
<b>Magnetic Agarose</b>	IP	About 50 $\mu\text{m}$	<p>Agarose beads Magnetic particles</p>		<ul style="list-style-type: none"> <li>-Large amount of IgG/per beads</li> <li>-Less sample loss because of the magnet</li> <li>-Easy to see the gel</li> <li>-Magnetic rack will be needed</li> </ul>	<ul style="list-style-type: none"> <li>★ Consider the volume of yield</li> <li>★ Want an easier procedure</li> </ul>
<b>Magnetic beads</b>	IP Purification	About 1.5 $\mu\text{m}$			<ul style="list-style-type: none"> <li>-Small amount of IgG/per beads</li> <li>-Almost no loss because of the magnet</li> <li>-Easy to disperse and can be used for screening</li> <li>-Easy to see the beads</li> <li>-Magnetic rack will be needed</li> </ul>	<ul style="list-style-type: none"> <li>★ Need smaller elution volume</li> <li>★ Want an easier procedure</li> </ul>
<b>Agarose</b>	IP Purification	About 100 $\mu\text{m}$			<ul style="list-style-type: none"> <li>-Large amount of IgG/per beads</li> <li>-Might lose during the wash</li> <li>-Hard to see the gel</li> <li>-No need of magnetic rack</li> </ul>	<ul style="list-style-type: none"> <li>★ Consider the volume of yield</li> <li>★ Want to reduce the cost</li> </ul>

### Example for the **Smart-IP** by Immunoprecipitation


#### ■ Anti-DDDDK-tag mAb (M185-10)



Lane 1: 293T cell lysate  
Lane 2: N-terminal DDDDK-tagged protein X (1  $\mu\text{g}$ ) in 293T lysate  
Lane 3: Internal DDDDK-tagged GFP (1  $\mu\text{g}$ ) in 293T lysate  
Lane 4: C-terminal DDDDK-tagged protein X (1  $\mu\text{g}$ ) in 293T lysate

Immunoblotted with anti-DDDDK-tag mAb-HRP-Direct (M185-7)


#### ■ Anti-HA-tag mAb (M180-10)



Lane 1: 293T  
Lane 2: HA-tagged protein expressed in 293T

Immunoblotted with M180-7


#### ■ Anti-HA-tag mAb (M132-11)



Sample: HA-tagged  $\beta$ -galactosidase/293T  
Lane 1: Input (10  $\mu\text{L}$ /lane)  
Lane 2: Post-IP beads of Anti-HA-tag mAb (M132-11)




■ Anti-His-tag mAb (D291-11)




Sample: HA-tagged β-galactosidase/293T  
Lane 1: Input (5 μL/lane)  
Lane 2: Post-IP beads of Anti-HA-tag mAb (D291-11)

■ Anti-V5-tag mAb (M167-10)




Lane 1: Insect medium  
Lane 2: V5-tagged TPO in insect medium  
Immunoblotted with PM003-7

■ Anti-RFP mAb (M165-11)




Sample: 293T cell lysate from 3x10<sup>9</sup> cells + mCherry protein\* 10 μg  
Lane 1: Input (10 μL/lane)  
Lane 2: Post-IP beads of Anti-RFP mAb (D291-11)  
\*Sample was provided by RIKEN.

■ Anti-Myc-tag mAb (M047-10)




Lane 1: 293T  
Lane 2: Myc-tagged protein expressed in 293T  
Immunoblotted with M192-7

■ Anti-GFP mAb (D153-10)



Lane 1: 293T  
Lane 2: GFP-fusion protein expressed in 293T  
Immunoblotted with 598-7

■ Anti-E-tag mAb (M198-9)



Lane 1: IP with Mouse IgG2a-Magnetic beads (M076-11)  
Lane 2: IP with Anti-E-tag mAb-Magnetic beads (M198-9)  
Immunoblotted with PM070

Smart-IP Tag antibody-conjugated magnetic beads or magnetic agarose

Code No.	Description	Clone	Isotype	Applications	Size
<b>Magnetic Agarose</b>					
M185-10	Anti-DDDDK-tag mAb-Magnetic Agarose	FLA-1	Mouse IgG2a <sub>κ</sub>	IP	100 tests (Gel: 1 mL)
M132-10	Anti-HA-tag mAb-Magnetic Agarose	5D8	Mouse IgG1 <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
M180-10	Anti-HA-tag mAb-Magnetic Agarose	TANA2	Mouse IgG2b <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
D291-10	Anti-His-tag mAb-Magnetic Agarose	OGHis	Mouse IgG2a <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
M047-10	Anti-Myc-tag mAb-Magnetic Agarose	PL14	Mouse IgG1 <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
M167-10	Anti-V5-tag mAb-Magnetic Agarose	1H6	Mouse IgG2a <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
D153-10	Anti-GFP mAb-Magnetic Agarose	RQ2	Rat IgG2a <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
M198-10	Anti-E-tag mAb-Magnetic Agarose	21D11	Mouse IgG2a <sub>κ</sub>	IP	20 tests (Gel: 200 μL)
<b>Magnetic Beads</b>					
M185-11	Anti-DDDDK-tag mAb-Magnetic Beads	FLA-1	Mouse IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M132-11	Anti-HA-tag mAb-Magnetic Beads	5D8	Mouse IgG1 <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M180-11	Anti-HA-tag mAb-Magnetic Beads	TANA2	Mouse IgG2b <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
D291-11	Anti-His-tag mAb-Magnetic Beads	OGHis	Mouse IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M047-11	Anti-Myc-tag mAb-Magnetic Beads	PL14	Mouse IgG1 <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M167-11	Anti-V5-tag mAb-Magnetic Beads	1H6	Mouse IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M215-11	Anti-V5-tag mAb-Magnetic Beads	OZA3	Mouse IgG2b <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
D153-11	Anti-GFP (Green Fluorescent Protein) mAb-Magnetic Beads	RQ2	Rat IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M165-11	Anti-RFP mAb-Magnetic Beads	3G5	Mouse IgG1 <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M198-9	Anti-E-tag mAb-Magnetic Beads	21D11	Mouse IgG2a <sub>κ</sub>	IP	20 tests (Slurry: 1 mL)
M075-11	Mouse IgG1 (isotype control)-Magnetic Beads	2E12	Mouse IgG1 <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M076-11	Mouse IgG2a (isotype control)-Magnetic Beads	6H3	Mouse IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M077-11	Mouse IgG2b (isotype control)-Magnetic Beads	3D12	Mouse IgG2b <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
M081-11	Rat IgG2a (isotype control)-Magnetic Beads	2H3	Rat IgG2a <sub>κ</sub>	IP, Purification	20 tests (Slurry: 1 mL)
MJS002V2	Protein G-Magnetic Beads	—	—	IP	10 mL (1% Slurry)
<b>Magnetic Rack</b>					
3190	Magnetic Rack				1.5 mL x 8 tubes

DDDDK-tag  
 HA-tag  
 His-tag  
 Myc-tag  
 V5-tag  
 mini-AID-tag  
 Fluorescent Protein Antibodies  
 Other Tag Antibodies  
 HRP-Direct  
 Smart-IP  
 Purification Kit and Gel  
 FAOs  
 Isotype Control Antibodies  
 Loading Control Antibodies  
 Organelle Marker

## Purification Kit & Gel

### PURIFICATION KIT

To purify target epitope-tagged proteins without any loss of protein activity and denature, all procedure during purification should be conducted under a neutral pH condition. Severe conditions such as acidic or alkaline elution can occasionally alter protein structures or biological functions of target proteins.


MBL's tagged Protein PURIFICATION KIT employs competitive polypeptides to elute target tagged proteins under a neutral pH for maintaining target protein function and conformation. Epitope tag specific antibodies in this kit are highly optimized for target recovery with high yield and purity from crude samples like cell lysates or culture supernatant of mammal transfectant cells.



\*This product is a set of purification gel, elution peptide, wash buffer, and mini column.

- ⊙ Recommended to who wants to purify from the small samples.
- ⊙ Almost everything is in this kit to start. Easy to get ready.
- ⊙ Perfect for pull-down assay.

#### Procedure Summary





### Purified in native conformations

Our purification procedure under a neutral pH condition enables maintaining target protein function and conformation. Useful/available for pull down assay as well to analyze co-immunoprecipitated proteins.


#### Purification and enzymatic activity of N-terminal His-tagged $\beta$ -galactosidase

His-tagged Protein PURIFICATION KIT (Code No. 3310)



#### Pull-down assay

HA-tagged Protein PURIFICATION KIT (Code No. 3320)




Interaction of p62/SQSTM1 with LC3 was confirmed by pull-down assay using the cell lysate of the HA-LC3 transfectant.

## High purity

Able to recover high purity target tagged proteins by using our original proprietary antibodies that have superior affinity and specificity to the epitope tag.

### ■ Purification of C-terminal His-tagged protein X from culture supernatant

His-tagged Protein PURIFICATION KIT (Code No. 3310)



## Wide variety of products

Wide range of purification kits for each tag-proteins.



## PURIFICATION GEL with Elution Peptide




\*This product is a set of purification gel and elution peptide.

◎ Recommended for the large scale of protein purification.

### ■ Purification of C-terminal His-tagged protein X from culture supernatant

His-tagged Protein PURIFICATION GEL with Elution Peptide (Code No. 3311)



SDS-PAGE (Coomassie Brilliant Blue Staining)  
Stable transfectant of CHO (Chinese Hamster Ovary) cells expressing C-terminal His-tagged protein were cultured for 7 days in DMEM medium containing 10% Fetal bovine serum. C-terminal His-tagged protein was purified from 50 mL of cultured medium according to the preceding protocol. Column bed volume was 0.5 mL. Elution was carried out with 2.5 mL of 1 mg/mL His tag peptide. Each fraction was 0.5 mL.

DDDDK-tag	
HA-tag	
His-tag	
Myc-tag	
V5-tag	
mini-AID-tag	
Fluorescent Protein Antibodies	
Other Tag Antibodies	
HRP-Direct	
Smart-IP	
Purification Kit and Gel	
FAOs	
Isotype Control Antibodies	
Loading Control Antibodies	
Organelle Marker Antibodies	

## PURIFICATION GEL




\*This product consists of a purification gel only.

- ⊙ Recommended for the large scale of protein purification.
- ⊙ For who wants to prepare elution peptide by themselves.

### ■ Purification of N-terminal DDDDK-tagged $\beta$ -galactosidase


DDDDK-tagged Protein PURIFICATION GEL (Code No. 3328R)



SDS-PAGE (Coomassie Brilliant Blue Staining)

Human embryonic kidney cell (293T) were transfected with pcDNA-DDDDK-tagged  $\beta$ -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (10 mL/100-mm dish x5) and purified according to the preceding protocol. Column bed volume was 0.25 mL. Elution was carried out with 2 mL of 0.1 mg/mL DDDDK-tag peptide. Each fraction was 0.25 mL.

## PURIFICATION CARTRIDGE





- ⊙ Ready-to-use-cartridge for the isolation of c-Myc-tagged protein
- ⊙ The cartridge can be used directly on automated ÄKTA or FPLC systems. (Cartridge is 1/16 inch female end)\*

\* Maximum working pressure: 0.3 MPa  
 \* ÄKTA is a trade mark of GE Healthcare (Life Sciences).

### ■ Example of Purification Result

c-Myc-tagged Protein PURIFICATION CARTRIDGE (Code No. 3306K)



Column bed volume : 1 mL  
 Sample : c-Myc-tagged GST/293T ( $5 \times 10^7$  cells)  
 Cell lysis buffer : 10 mM Tris-HCl, 150 mM NaCl, 1% NP-40 (pH 7.5)  
 Wash buffer : 0.05% Tween/PBS (pH 7.2)  
 Elution buffer : 0.5 mg/ml c-Myc peptide in PBS  
 Flow rate : 0.75 mL/min  
 Chromatography system : ÄKTA explorer 10S (GE Healthcare)


## Tagged Protein Magnetic Purification Kits



\*Sold separately Code No. 3190


- Used magnet for our popular purification kits.
- Magnetic conjugated antibodies, elution peptide, wash solutions are in package.
- Recommended for purifying small amount samples.
- Set of needed reagents.
- No need to centrifuge.
- High yield
- Maintaining target protein function and conformation

### ■ Procedure Summary



### ■ c-Myc-tagged GST protein purification


c-Myc-tagged Protein Magnetic Purification Kit (Code No. 3340)



Lane1: Input  
 Lane2: Eluate  
 Sample: c-Myc-tagged GST/*E.coli* JM109

### ■ V5-tagged $\beta$ -galactosidase purification

V5-tagged Protein Magnetic Purification Kit (Code No. 3341)



Lane 1: Input  
 Lane 2: Eluate  
 Sample: V5-tagged  $\beta$ -galactosidase/HEK293T

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAQs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies



Try out our trial kits!

Code No.	Description	Size
3325	DDDDK-tagged Protein PURIFICATION KIT	20 tests
3325A	DDDDK-tagged Protein PURIFICATION KIT (Trial Kit)	2 tests
3326	DDDDK-tagged Protein PURIFICATION GEL with Elution Peptide	gel 1 mL, peptide 5 mg
3327	DDDDK-tagged Protein PURIFICATION GEL with Elution Peptide	gel 5 mL, peptide 25 mg
3328R	DDDDK-tagged Protein PURIFICATION GEL	gel 5 mL
3329R	DDDDK-tagged Protein PURIFICATION GEL	gel 25 mL
3326KR	DDDDK-tagged Protein PURIFICATION CARTRIDGE	1 mL x 1
3325-205	DDDDK-tag peptide (DYKDDDDK)	1 mg x 5
3343	DDDDK-tagged Protein Magnetic Purification Kit	1 kit
3343A	DDDDK-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3320	HA-tagged Protein PURIFICATION KIT	20 tests
3320A	HA-tagged Protein PURIFICATION KIT Trial Kit	2 tests
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide (YPYDVPDYA)	2 mg x 5
3342	HA-tagged Protein Magnetic Purification Kit	1 kit
3342A	HA-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3310	His-tagged Protein PURIFICATION KIT	20 tests
3310A	His-tagged Protein PURIFICATION KIT Trial Kit	2 tests
3311	His-tagged Protein PURIFICATION GEL	gel 1 mL x 1, peptide 2 mg x 5
3312	His-tagged Protein PURIFICATION GEL	gel 1 mL x 5, peptide 2 mg x 25
3310-205	His-tag peptide (XXX-(6xHis)-XXX)	2 mg x 5
3305	c-Myc-tagged Protein MILD PURIFICATION KIT Ver.2	20 tests
3305A	c-Myc-tagged Protein MILD PURIFICATION KIT Ver.2 (Trial Kit)	2 tests
3306	c-Myc-tagged Protein MILD PURIFICATION GEL	gel 1 mL, peptide 1 mg
3307	c-Myc-tagged Protein MILD PURIFICATION GEL	gel 1 mL x 5, peptide 1 mg x 5
3306K	c-Myc-tagged Protein PURIFICATION CARTRIDGE	1 mL x 1
3300-205	c-Myc-tag peptide (EQKLISEEDL)	1 mg x 5
3340	c-Myc-tagged Protein Magnetic Purification Kit	1 kit
3340A	c-Myc-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3317	V5-tagged Protein Purification Kit Ver.2	20 tests
3317A	V5-tagged Protein Purification Kit Ver.2 (Trial Kit)	2 tests
3318	V5-tagged Protein Purification Gel Ver.2	1 mL
3315-205	V5-tag peptide (GKPIPPELLGLDST)	2 mg x 5
3341	V5-tagged Protein Magnetic Purification Kit	1 kit
3341A	V5-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3190	Magnetic Rack	1.5 mL x 8 tubes

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAQs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

### 1. I cannot purify epitope-tagged target proteins from my samples.

- If you can successfully confirm that the target protein was expressed rightly by other methods (e.g. western blotting), check whether the epitope-tagged target protein binds to anti-epitope tag antibody conjugated beads (tag beads) or not.
- Briefly, after incubation the tag beads with the samples, please add SDS-PAGE sample buffer to tag beads, boil it for 5 minutes and perform SDS-PAGE or western blotting.
- If you cannot observe the target band as the results, please read the Q2.
- If you can obtain the target band as the results, please read the Q3.

### 2. My epitope-tagged target protein does not bind to the tag beads.

- There are several possible causes as follows;
  - Your target protein is insoluble or aggregated.
  - Your lysis buffer contains guanidine and/or high concentrated urea. Some reagents can inhibit antigen-antibody reaction.
- Please refer to “Additional information” in the data sheet for each PURIFICATION KIT. “Additional information” is about the propriety of usage of reagents contained in cell lysis buffer.
- Sometimes the compatibility between epitope-tags and target proteins affects to protein solubility. In such cases, try to change the location of epitope tag or use other epitope tags.

### 3. My epitope-tagged target protein binds to tag beads, but I cannot elute the protein with peptide.

- Please try to prolong the incubation time or increase elution peptide volume.
- If your proteins are easy to be aggregated, please change buffer composition.
- e.g. Use lysate buffer for washing step instead of Wash Solution which is included in the PURIFICATION KIT.
- When you change composition of buffers, please refer to “Additional information” in the data sheet for each PURIFICATION KIT. “Additional information” is about the propriety of usage of reagents contained in cell lysis buffer.

### 4. Which lysis buffer should I use?

- The buffer which contains usual detergent (e.g. NP-40 and Tween 20) and 0.1%-SDS can be used.

### 5. Does the epitope tag location affect purification efficiency?

- In many cases, MBL's PURIFICATION KITS can purify the epitope-tagged target proteins regardless of tag location. However, preliminary study is important in any cases.

### 6. Can I elute the epitope-tagged target protein at 4°C ?

- Most of MBL's PURIFICATION KITS accept 4°C in elution step. However, in His tagged Protein PURIFICATION KIT, the purification efficiency is reduced to almost halve if you elute the epitope-tagged target protein at 4°C . If you have to perform elution step at 4°C , please incubate His-tag beads with Elution Peptide Solution at 4°C overnight before elution.

DDDDK-tag	HA-tag	His-tag	Myc-tag	V5-tag	mini-AID-tag	Fluorescent Protein Antibodies	Other Tag Antibodies	HRP-Direct	Smart-IP	Purification Kit and Gel	FAQs	Isotype Control Antibodies	Loading Control Antibodies	Organelle Marker Antibodies
-----------	--------	---------	---------	--------	--------------	--------------------------------	----------------------	------------	----------	--------------------------	------	----------------------------	----------------------------	-----------------------------

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
<b>FAQs</b>
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies

## 7. Can I purify the epitope-tagged target protein from inclusion bodies of *E. coli* ?


→ Some buffers which contain guanidine and/or high concentrated urea cannot be used for purification using MBL's PURIFICATION KITS. Aggregated or insoluble proteins cannot be purified, because such proteins cannot bind to tag beads. Therefore, please establish appropriate buffer condition with referring to "Additional information" in the data sheet for each PURIFICATION KIT. "Additional information" is about the propriety of usage of reagents contained in cell lysis buffer.

## 8. Can I elute the epitope-tagged target protein using buffers other than Elution Peptide Solution?

→ Yes you can. SDS-PAGE sample buffer, acid elution solution, and alkali elution solution can be used. However, the target proteins may lose protein activity and native conformation if you use such elution buffers.

Acid elution solution: 0.1 M Glycine-HCl, pH 3.0 (Neutralize the elution immediately with 1 M Tris-HCl, pH 8.0)

Alkali elution solution: 0.1M NH<sub>3</sub>, pH 11.3 (Neutralize the elution immediately with 1 N acetic acid)





## Related Antibodies

### Isotype Control Antibodies – Negative Control Antibodies

Code No.	Description	Clone	Isotype	Applications	Size
M075-3	Mouse IgG1	2E12	Mouse IgG1 $\kappa$	IP, FCM	100 $\mu$ g/100 $\mu$ L
M075-3M2	Mouse IgG1 (Functional Grade)	2E12	Mouse IgG1 $\kappa$	–	500 $\mu$ g
M075-4	Mouse IgG1-FITC	2E12	Mouse IgG1 $\kappa$	FCM	50 $\mu$ g/1 mL
M075-5	Mouse IgG1-PE	2E12	Mouse IgG1 $\kappa$	FCM	1 mL (50 tests)
M075-6	Mouse IgG1-Biotin	2E12	Mouse IgG1 $\kappa$	FCM	50 $\mu$ g/50 $\mu$ L
M075-8	Mouse IgG1-Agarose <span style="background-color: #0070C0; color: white; padding: 2px;">Agarose</span>	2E12	Mouse IgG1 $\kappa$	IP	400 $\mu$ L
M075-11	Mouse IgG1-Magnetic Beads <span style="background-color: #0070C0; color: white; padding: 2px;">Smart-IP</span>	2E12	Mouse IgG1 $\kappa$	IP	20 tests (Slurry: 1 mL)
M075-12	Mouse IgG1 (isotype control)-ALP	2E12	Mouse IgG1 $\kappa$	ELISA	50 $\mu$ L
M075-A48	Mouse IgG1-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	2E12	Mouse IgG1 $\kappa$	FCM	100 $\mu$ g
M075-A64	Mouse IgG1-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	2E12	Mouse IgG1 $\kappa$	FCM	100 $\mu$ g
M076-3	Mouse IgG2a	6H3	Mouse IgG2a $\kappa$	IP, FCM	100 $\mu$ g/100 $\mu$ L
M076-3M2	Mouse IgG2a (Functional Grade)	6H3	Mouse IgG2a $\kappa$	–	500 $\mu$ g/100 $\mu$ L
M076-4	Mouse IgG2a-FITC	6H3	Mouse IgG2a $\kappa$	FCM	1 mL
M076-5	Mouse IgG2a-PE	6H3	Mouse IgG2a $\kappa$	FCM	1 mL (50 tests)
M076-6	Mouse IgG2a-Biotin	6H3	Mouse IgG2a $\kappa$	FCM	50 $\mu$ g/50 $\mu$ L
M076-11	Mouse IgG2a-Magnetic Beads <span style="background-color: #0070C0; color: white; padding: 2px;">Smart-IP</span>	6H3	Mouse IgG2a $\kappa$	IP	20 tests (Slurry: 1 mL)
M076-12	Mouse IgG2a (isotype control)-ALP	6H3	Mouse IgG2a $\kappa$	ELISA	50 $\mu$ L
M076-A48	Mouse IgG2a-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	6H3	Mouse IgG2a $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M076-A64	Mouse IgG2a-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	6H3	Mouse IgG2a $\kappa$	FCM	100 $\mu$ L/1 mg/mL
M077-3	Mouse IgG2b	3D12	Mouse IgG2b $\kappa$	IP, FCM	100 $\mu$ g/100 mL
M077-3M2	Mouse IgG2b (Functional Grade)	3D12	Mouse IgG2b $\kappa$	–	500 $\mu$ g/100 $\mu$ L
M077-4	Mouse IgG2b-FITC	3D12	Mouse IgG2b $\kappa$	FCM	50 $\mu$ g/1 mL
M077-5	Mouse IgG2b-PE	3D12	Mouse IgG2b $\kappa$	FCM	1 mL (50 tests)
M077-6	Mouse IgG2b-Biotin	3D12	Mouse IgG2b $\kappa$	FCM	50 $\mu$ g/50 $\mu$ L
M077-11	Mouse IgG2b-Magnetic Beads <span style="background-color: #0070C0; color: white; padding: 2px;">Smart-IP</span>	3D12	Mouse IgG2b $\kappa$	IP	20 tests (Slurry: 1 mL)
M077-12	Mouse IgG2b (isotype control)-ALP	3D12	Mouse IgG2b $\kappa$	ELISA	50 $\mu$ L
M077-A48	Mouse IgG2b-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	3D12	Mouse IgG2b $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M077-A64	Mouse IgG2b-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	3D12	Mouse IgG2b $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M078-3	Mouse IgG3	6A3	Mouse IgG3	IP, FCM	100 $\mu$ g/100 $\mu$ L
M078-3M2	Mouse IgG3 (Functional Grade)	6A3	Mouse IgG3	–	500 $\mu$ g/100 $\mu$ L
M078-4	Mouse IgG3-FITC	6A3	Mouse IgG3	FCM	50 $\mu$ g
M078-6	Mouse IgG3 -Biotin	6A3	Mouse IgG3	FCM	50 $\mu$ g/50 $\mu$ L
M079-3	Mouse IgM	7E10	Mouse IgM	FCM	100 $\mu$ g/100 $\mu$ L
M080-3	Rat IgG1	1H5	Rat IgG1	IP, FCM	100 $\mu$ g/100 $\mu$ L
M080-3M2	Rat IgG1 (Functional Grade)	1H5	Rat IgG1	–	500 $\mu$ g/100 $\mu$ L
M080-4	Rat IgG1-FITC	1H5	Rat IgG1	FCM	50 $\mu$ g
M080-5	Rat IgG1-PE	1H5	Rat IgG1	FCM	1 mL (50 tests)
M080-A48	Rat IgG1-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	1H5	Rat IgG1	FCM	100 $\mu$ g
M080-A64	Rat IgG1-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	1H5	Rat IgG1	FCM	100 $\mu$ g/100 $\mu$ L
M081-3	Rat IgG2a	2H3	Rat IgG2a $\kappa$	IP, FCM	100 $\mu$ g/100 $\mu$ L
M081-3M2	Rat IgG2 (Functional Grade)	2H3	Rat IgG2a $\kappa$	–	500 $\mu$ g/100 $\mu$ L
M081-4	Rat IgG2a-FITC	2H3	Rat IgG2a $\kappa$	FCM	50 $\mu$ g/1 mL
M081-5	Rat IgG2a-PE	2H3	Rat IgG2a $\kappa$	FCM	1 mL (50 tests)
M081-8	Rat IgG2a-Agarose <span style="background-color: #0070C0; color: white; padding: 2px;">Agarose</span>	2H3	Rat IgG2a $\kappa$	IP	20 tests (Gel: 200 $\mu$ L)
M081-11	Rat IgG2a-Magnetic Beads <span style="background-color: #0070C0; color: white; padding: 2px;">Smart-IP</span>	2H3	Rat IgG2a $\kappa$	IP	20 tests (Slurry: 1 mL)
M081-A48	Rat IgG2a-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	2H3	Rat IgG2a $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M081-A64	Rat IgG2a-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	2H3	Rat IgG2a $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M082-3	Rat IgG2c	6E12	Rat IgG2c	IP, FCM	100 $\mu$ g/100 $\mu$ L
M082-3M2	Rat IgG2a (Functional Grade)	6E12	Rat IgG2c	–	500 $\mu$ g/100 $\mu$ L
M082-4	Rat IgG2c-FITC	6E12	Rat IgG2c	FCM	50 $\mu$ g
M090-3	Rat IgG2b	3G8	Rat IgG2b $\kappa$	IP, FCM	100 $\mu$ g
M090-3M2	Rat IgG2b (Functional Grade)	3G8	Rat IgG2b $\kappa$	–	500 $\mu$ g
M090-4	Rat IgG2b-FITC	3G8	Rat IgG2b $\kappa$	FCM	1 mL
M090-5	Rat IgG2b-PE	3G8	Rat IgG2b $\kappa$	FCM	1 mL (50 tests)
M090-A48	Rat IgG2b-Alexa Fluor <sup>®</sup> 488 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	3G8	Rat IgG2b $\kappa$	FCM	100 $\mu$ g
M090-A64	Rat IgG2b-Alexa Fluor <sup>®</sup> 647 <span style="background-color: #90EE90; border: 1px solid black; padding: 2px;">Alexa Fluor<sup>®</sup></span>	3G8	Rat IgG2b $\kappa$	FCM	100 $\mu$ g/100 $\mu$ L
M189-3	Syrian Hamster IgG	ttko	Hamster IgG	IP, FCM	100 $\mu$ g/100 $\mu$ L
M199-3	Armenian Hamster IgG	ttko2	Hamster IgG	IP, FCM	100 $\mu$ g/100 $\mu$ L
PM035	Normal Rabbit IgG	Polyclonal	Rabbit IgG	IP, FCM	500 $\mu$ g/100 $\mu$ L
PM035-8	Normal Rabbit IgG-Agarose <span style="background-color: #0070C0; color: white; padding: 2px;">Agarose</span>	Polyclonal	Rabbit IgG	IP	Gel: 200 $\mu$ L
PM067	Normal Guinea Pig IgG	Polyclonal	Guinea Pig IgG	IP, FCM	100 $\mu$ g/100 $\mu$ L
PM084	Normal Chicken IgY	Polyclonal	Chicken IgY	FCM	200 $\mu$ g/200 $\mu$ L
PM084-4	Normal Chicken IgY-FITC	Polyclonal	Chicken IgY	FCM	100 $\mu$ g/100 $\mu$ L
PM084-5	Normal Chicken IgY-PE	Polyclonal	Chicken IgY	FCM	1 mL (50 tests)
PM094	Normal Goat IgG	Polyclonal	Goat IgG	IC	500 $\mu$ g/100 $\mu$ L
M194-3	Human IgG1 isotype control chimeric mAb	2E12G1-2	Human IgG1	–	100 $\mu$ g/100 $\mu$ L
M195-3	Human IgG2 isotype control chimeric mAb	2E12G2-18	Human IgG2	–	100 $\mu$ g/100 $\mu$ L

DDDDK-tag
HA-tag
His-tag
Myc-tag
V5-tag
mini-AID-tag
Fluorescent Protein Antibodies
Other Tag Antibodies
HRP-Direct
Smart-IP
Purification Kit and Gel
FAOs
Isotype Control Antibodies
Loading Control Antibodies
Organelle Marker Antibodies


## Loading Control Antibodies

Code No.	Description	molecular weight	Clone	Isotype	Applications	Size
M171-3	Anti-GAPDH mAb	around 35 kDa (36 kDa)	3H12	Mouse IgG2a $\kappa$	WB	300 $\mu$ g/100 $\mu$ L
M171-7	Anti-GAPDH mAb-HRP-Direct <b>HRP-Direct</b>	around 35 kDa (36 kDa)	3H12	Mouse IgG2a $\kappa$	WB	50 $\mu$ L
M177-3	Anti- $\beta$ -Actin mAb	around 40 kDa (42 kDa)	6D1	Mouse IgG1 $\kappa$	WB, IP	100 $\mu$ g/100 $\mu$ L
PM053	Anti- $\beta$ -Actin pAb	around 40 kDa (42 kDa)	Polyclonal	Rabbit Ig (aff.)	WB, IP	100 $\mu$ L
PM053-7	Anti- $\beta$ -Actin pAb-HRP-Direct <b>HRP-Direct</b>	around 40 kDa (42 kDa)	Polyclonal	Rabbit Ig (aff.)	WB	50 $\mu$ L
M175-3	Anti- $\alpha$ -Tubulin mAb	around 50 kDa (50 kDa)	2F9	Mouse IgG2a $\kappa$	WB, IP, IC	200 $\mu$ g/100 $\mu$ L
PM054	Anti- $\alpha$ -Tubulin pAb	around 50 kDa (50 kDa)	Polyclonal	Rabbit Ig (aff.)	WB, IP, IC	100 $\mu$ L
PM054-7	Anti- $\alpha$ -Tubulin pAb-HRP-Direct <b>HRP-Direct</b>	around 50 kDa (50 kDa)	Polyclonal	Rabbit Ig (aff.)	WB	50 $\mu$ L
PM064	Anti-Lamin B1 pAb	around 70 kDa (66 kDa)	Polyclonal	Rabbit Ig (aff.)	WB, IP, IC	100 $\mu$ L
PM088	Anti-Vinculin pAb	around 120 kDa (124 kDa)	Polyclonal	Rabbit Ig (aff.)	WB	100 $\mu$ L
8480	Loading Control Set (HRP-Direct) (Contains M171-7, PM054-7, PM053-7) <b>HRP-Direct</b>				WB	3 vials (50 $\mu$ L each)

## Organelle Marker Antibodies

Code No.	Description	Clone	Isotype	Target	Applications	Size
PM062	Anti-EEA1 pAb	Polyclonal	Rabbit Ig (aff.)	Early endosome	WB, IP, IC	100 $\mu$ L
M176-3	Anti-EEA1 mAb	3C10	Mouse IgG2a $\kappa$	Early endosome	WB, IP, IC	100 $\mu$ g/100 $\mu$ L
M176-A48	Anti-EEA1 mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	3C10	Mouse IgG2a $\kappa$	Early endosome	IC	100 $\mu$ g/100 $\mu$ L
M176-A59	Anti-EEA1 mAb-Alexa Fluor <sup>®</sup> 594 <b>Alexa Fluor<sup>®</sup></b>	3C10	Mouse IgG2a $\kappa$	Early endosome	IC	100 $\mu$ g/100 $\mu$ L
M176-A64	Anti-EEA1 mAb-Alexa Fluor <sup>®</sup> 647 <b>Alexa Fluor<sup>®</sup></b>	3C10	Mouse IgG2a $\kappa$	Early endosome	IC	100 $\mu$ g/100 $\mu$ L
PM060	Anti-Calnexin pAb	Polyclonal	Rabbit Ig (aff.)	Endoplasmic reticulum	WB, IP, IC	100 $\mu$ L
M178-3	Anti-Calnexin mAb	4F10	Mouse IgG2a $\kappa$	Endoplasmic reticulum	WB, IP, IC	100 $\mu$ g/100 $\mu$ L
PM059	Anti-KDEL pAb	Polyclonal	Rabbit Ig (aff.)	Endoplasmic reticulum	WB, IC, IH*	100 $\mu$ L
M181-3	Anti-KDEL mAb	1D5	Mouse IgG2a $\kappa$	Endoplasmic reticulum	WB, IP*, IC, IH*	100 $\mu$ g/100 $\mu$ L
PM061	Anti-GM130 pAb	Polyclonal	Rabbit Ig (aff.)	Golgi body	WB, IP, IC	100 $\mu$ L
M179-3	Anti-GM130 mAb	5G8	Mouse IgG2a $\kappa$	Golgi body	WB, IP, IC	100 $\mu$ g/100 $\mu$ L
M179-A48	Anti-GM130 mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	5G8	Mouse IgG2a $\kappa$	Golgi body	IC	100 $\mu$ g/100 $\mu$ L
PM063	Anti-COX4 pAb	Polyclonal	Rabbit Ig (aff.)	Mitochondria	WB, IP, IC	100 $\mu$ L
PM064	Anti-Lamin B1 pAb	Polyclonal	Rabbit Ig (aff.)	Nuclear membrane	WB, IP, IC	100 $\mu$ L
D115-3	Anti-CENP-A mAb	3-19	Mouse IgG1	Centromere	WB, IC, IH, ChIP*	100 $\mu$ g/100 $\mu$ L
PD030	Anti-CENP-C pAb	Polyclonal	Guinea pig IgG	Centromere	WB, IP, IC	100 $\mu$ L
PM054	Anti- $\alpha$ -Tubulin pAb	Polyclonal	Rabbit Ig (aff.)	Microtubule	WB, IP, IC	100 $\mu$ L
PM054-7	Anti- $\alpha$ -Tubulin pAb-HRP-Direct <b>HRP-Direct</b>	Polyclonal	Rabbit Ig (aff.)	Microtubule	WB	50 $\mu$ L
M175-3	Anti- $\alpha$ -Tubulin mAb	2F9	Mouse IgG2a $\kappa$	Microtubule	WB, IP, IC	200 $\mu$ g/100 $\mu$ L
M175-A48	Anti- $\alpha$ -Tubulin mAb-Alexa Fluor <sup>®</sup> 488 <b>Alexa Fluor<sup>®</sup></b>	2F9	Mouse IgG2a $\kappa$	Microtubule	IC	100 $\mu$ g/100 $\mu$ L
PD033	Anti- $\beta$ 1-Tubulin pAb	Polyclonal	Rabbit Ig (aff.)	Microtubule	WB, IC	100 $\mu$ L
PM036	Anti-LC3 pAb	Polyclonal	Rabbit IgG	Autophagosome	WB, IP, FCM, IC, IH	100 $\mu$ L
PD014	Anti-LC3 pAb	Polyclonal	Rabbit IgG	Autophagosome	WB, IC*, IH*	100 $\mu$ L
M152-3	Anti-LC3 mAb	4E12	Mouse IgG1 $\kappa$	Autophagosome	WB, IP, FCM, IC, IH*, Immuno-EM	200 $\mu$ g/100 $\mu$ L
M186-3	Anti-LC3 mAb	8E10	Mouse IgG2a $\kappa$	Autophagosome	WB	100 $\mu$ g/100 $\mu$ L
M186-7	Anti-LC3 mAb-HRP-Direct <b>HRP-Direct</b>	8E10	Mouse IgG2a $\kappa$	Autophagosome	WB	50 $\mu$ L

## Organelle Marker Antibodies



DDDDK-tag	
HA-tag	
His-tag	
Myc-tag	
V5-tag	
mini-AID-tag	
Fluorescent Protein Antibodies	
Other Tag Antibodies	
HRP-Direct	
Smart-IP	
Purification Kit and Gel	
FAQs	
Isotype Control Antibodies	
Loading Control Antibodies	
<b>Organelle Marker Antibodies</b>	

Produced by

**MBL** MEDICAL & BIOLOGICAL  
LABORATORIES CO., LTD.

A JSR Life Sciences Company

SUMITOMO FUDOSAN SHIBADAIMON NICHOME BLDG.  
2-11-8 Shibadaimon, Minato-ku, Tokyo 105-0012 Japan  
TEL: +81-3-6854-3614 E-mail: support@mbi.co.jp  
URL: <https://www.mblbio.com/bio/g/>